


Lactosucrose has a sweet taste similar to sugar and regulates the functions of the intestines. In Japan, food products that contain appropriate amount of lactosucrose can be approved and labeled as foods for specified health uses (FOSHU). In the standard regulation system of FOSHU, HPLC method using silica gel chemically bonded with carbamoyl groups is described as the determination method for lactosucrose.

In this note, Inertsil Amide was used as a separation column for the analysis. Concentration of lactosucrose in commercially available beverage was successfully determined.

(K. Kanno)

A Chromatogram Obtained from Standard Solution



HPLC conditions

Column : Inertsil Amide
(5 μ m, 250 x 4.6 mm I.D.)

Eluent : A) CH₃CN
B) H₂O
A/B = 71/29, v/v

Flow rate : 1.2 mL/min


Col. Temp. : 35 °C

Detection : RI

Inj. Vol. : 20 μ L

* Aqueous solution of lactosucrose (10 mg/mL) was 2-fold diluted with acetonitrile.

Chemical Structure



1. Lactosucrose

Structures are created using Chemistry 4-D Draw which is provided by ChemInnovation Software, Inc.






Figure 1. Calibration curve for lactosucrose

An Example of Pretreatment (for liquid sample)





A Chromatogram Obtained from a Beverage Containing Lactosucrose



A know-how for using amide column ③

(A) water 100 % (B) water/acetonitrile = 50/50

When using amide column, it is preferable that sample is dissolved in acetonitrile. In general, peak shape is improved as acetonitrile concentration of sample solvent is increased. Chromatograms shown right are obtained from lactosucrose solution dissolved in water 100% (A) or water/acetonitrile = 50/50 (B).



1

Column : Inersil Amide
(5 μm, 250 x 4.6 mm I.D.)
Eluent : A) CH₃CN
B) H₂O
A/B = 71/29, v/v
Flow rate : 1.2 mL/min
Col. Temp. : 35 °C
Detection : RI
Inj. Vol. : 20 μL

GL Sciences disclaims any and all responsibility for any injury or damage which may be caused by this data directly or indirectly. We reserve the right to amend this information or data at any time and without any prior announcement.

GL Sciences, Inc. Japan

22-1 Nishishinjuku 6-Chome
Shinjuku-ku, Tokyo,
163-1130, Japan
Phone: +81-3-5323-6620
Fax: +81-3-5323-6621
Email: world@glsc.co.jp
Web: www.glsciences.com

GL Sciences B.V.

De Sleutel 9
5652 AS Eindhoven
The Netherlands
Phone: +31 (0)40 254 95 31
Email: info@glsciences.eu
Web: www.glsciences.eu

GL Sciences (ShangHai) Ltd.

Tower B, Room 2003,
Far East International Plaza,
NO.317 Xianxia Road,
Changning District.
Shanghai, China P.C. 200032
Phone: +86 (0)21-6278-2272
Email: contact@glsciences.com.cn
Web: www.glsciences.com.cn

GL Sciences, Inc. USA

4733 Torrance Blvd. Suite 255
Torrance, CA 90503
Phone: 310-265-4424
Fax: 310-265-4425
Email: info@glsciencesinc.com
Web: www.glsciencesinc.com

International Distributors

Visit our Website at:

<https://www.glsciences.com/company/distributor.html>