


Analysis of Febuxostat

(Under the Condition of the Japanese Pharmacopoeia 18th Supplement II, Febuxostat, Related substances (ii))


Test for required detectability

1-a. Febuxostat 50 μ g/L
 2-a. Related substances B* 50 μ g/L
 (Peak area of 1-a) / (Peak area of 1-b) : $(7 \leq) 7.34 (\leq 13)$
 (Peak area of 2-a) / (Peak area of 2-b) : $(7 \leq) 7.11 (\leq 13)$




System performance

1-b. Febuxostat 500 μ g/L
 2-b. Related substances B* 500 μ g/L
 Resolution (1.2) : $5.6 (\geq 2.0)$



System repeatability

1-c. Febuxostat 500 μ g/L
 RSD of the
 Peak area (%) (n=6) : $0.91 (\leq 2.0)$



Conditions

System : Chromaster HPLC system (HITACHI)
Column : InertSustain C30 (GL Sciences Inc.)
 (HP 3 μ m, 150 x 4.6 mm I.D.)
Column Cat. No. : 5020-17184
Eluent : A) 0.05% CF₃COOH in H₂O
 B) 0.05% CF₃COOH in CH₃CN
 A/B = 55/45, v/v
Flow Rate : 1.19 mL/min
Col. Temp. : 15 °C
Detection : UV 317 nm
Injection Vol. : 20 μ L
Sample : Standard

Analyte:

1. Febuxostat
 2. Related substances B*

*2-(4-Butoxy-3-cyanophenyl)-4-methyl-1,3-thiazole-5-carboxylic acid