


Sample Preparation

800 μL sample solution :
Add acetonitrile to PA
sugar chain sample
solution and adjust
the acetonitrile final
concentration from 90 to
95 %.

Centrifugation
Speed : 2,300 x g



1. Conditioning

Add 500 μL solution
mixed with 250 μL 0.1 %
formic acid* in water
and 250 μL 0.1 % formic
acid in acetonitrile
↓ Centrifuge
Add 500 μL solution
mixed with 50 μL 0.1 %
formic acid* in water
and 450 μL 0.1 % formic
acid in acetonitrile

2. Adsorption

Add 800 μL
sample solution

3. Rinsing

Add 500 μL
solution mixed
with 50 μL 0.1 %
formic acid* in
water and 450 μL
0.1 % formic acid
in acetonitrile


4. Elution

Add 50-800 μL
0.1 % formic acid in
50 % acetonitrile


**Purified
Sample**

* Acetic acid or TFA can
also be used as an
alternative to formic
acid.

Without MonoSpin NH2



Purification of PA using MonoSpin NH2



Conditions

Column : NH_2 Column (5 μm , 250 x 4.6 mm I.D.)
Eluent : A) H_2O /Acetonitrile = 5/95 0.1 % Formic acid
B) H_2O /Acetonitrile = 95/5 0.1 % Formic acid
A/B = 90/10-10 min-90/10-40 min-60/40
Flow Rate : 1 mL/min
Detection : FL Em 320 nm, Ex 400 nm
Injection Vol. : 1.5 μL

Purified PA sugar chain by HILIC mode.

MonoSpin NH2 additionally removes residual fluorescent labeling reagents.