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How to recover naturally sealed tissues in implantology ?

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
From its very beginnings and still today, oral implantology has been faced with the challenge of maintaining the integrity of peri-implant tissues, particularly those that make up the biological width.

It is well known that the tissues surrounding a natural tooth are different from those surrounding a dental implant. One particularity is the absence of the periodontal ligament between the alveolar bone socket and the implant, absence which interferes with and weakens tissue healing.

Good oral health means preserving all tissues, especially those attached to the tooth and to the implants. These attachments, although different in contact with a tooth or an implant, must nonetheless ensure the same shield function against inflammatory, infectious and occlusal stress.

Thanks to past research and recent publications, clinicians and scientists have gained a better understanding of the biological width and its physiology. New paradigms are starting to make consensus, such as implant bioforms, the choice of pure Titanium and Zirconia biomaterials, the machining accuracy of the implant connections, always with a view to protecting the integrity of the supporting tissues and to achieving a high long-term success rate of the implant treatment plans.

The challenge for both dental practitioners and implant manufacturers is to offer minimally invasive treatments, capable of taking into account all the factors influencing the protection and the stability of this biological width.

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

3.1 Implantologie Tissue Level Tissue Level implantology

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Accurate manufacturing tolerances of up to 5 µm are recommended for implant connections to ensure precision and stability of prosthetic assemblies (no micro-gap).

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The grade 4 Titanium body combines macro and micro roughness, providing biocompatibility and reduced healing times.

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
The Zirconia collar preserves the biological width by maintaining the attachment of connective fibers (horizontal and oblique) which support the tissues, inhibit invagination of the junctional epithelium and prevent gingival recession.

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The periodontal health thus preserved reflects the tight seal of the tissues around the implant :

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

Better results on 5 biological criteria


COMPARAISON DE L'INDICE GINGIVAL
COMPARISON OF GINGIVAL INDEX SCORE


COMPARAISON DE L'INDICE DE PLAQUE
COMPARISON OF PLAQUE INDEX SCORE


COMPARAISON DE LA PROFONDEUR DE SONDAGE PARODONTAL
COMPARISON OF PERIODONTAL POCKET DEPTH

COMPARAISON DE L'INDICE DE SAIGNEMENT
COMPARISON OF MODIFIED SULCULAR BLEEDING INDEX

COMPARAISON DE L'INDICE ESTHÉTIQUE
COMPARISON OF ESTHETIC VALUE (VAS)

S. Suryavanshi and M. Verma, "Comparative evaluation of the soft tissue response and esthetics of the titanium implants with Zirconium collar," Department of Prosthodontics, Crown and Bridge, Maulana Azad Institute of Dental Sciences, New Delhi, Clinical report, 2014.

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
By analyzing the ratio of Connective / Epithelial Tissue, surface-treated titanium represents the material of choice for complete osseointegration in a reduced time. As for the Zirconia surface, it offers a mucosal attachment similar to natural teeth.

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By combining these biomaterials, the resulting tight seal of the tissues to inflammatory and bacterial colonization enables to maintain the junctional epithelium at a distance from the bone, just like a natural tooth.

Sealed Tissue Recovery for Optimal Natural Growth

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In a study published in 2019, Korean and Japanese researchers not only analyzed the differential healing patterns of mucosal seal (in-vitro and in-vivo) on various biomaterials (Titanium, Zirconia, natural tooth), but also proposed a biologically more suitable type of implant.

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

To achieve naturally sealed tissues in implantology, their conclusion suggests :

✔ **Implant Tissue Level**
Tissue Level implant

✔ **Traitement de surface du corps en Titane**
Surface treatment of the Titanium body

✔ **Surmonté d'un col gingival en Zircone**
Featuring a gingival Zirconia collar

RATIO DE L'ESPACE BIOLOGIQUE
BIOLOGICAL WIDTH RATIO

D.-J. Lee, J.-S. Ryu, M. Shimono, K.-W. Lee, J.-M. Lee, and H.-S. Jung, "Differential Healing Patterns of Mucosal Seal on Zirconia and Titanium Implant," *Front. Physiol.*, vol. 10, no. 796, 2019, doi: 10.3389/fphys.2019.00796.

6 Pour quelles indications cliniques ? For which clinical indications ?

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
Cement-retained crowns and bridges

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


Custom-made prosthesis :
abutments and bars

7 Comment stimuler les tissus en phase de cicatrisation ? How to stimulate the tissues during the healing phase ?

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Healing crown

Without a periodontal ligament, the mechanical forces applied to the teeth are transferred to the alveolar bone, process which is influenced by the shape and material of the prosthetic restoration. In order to control these occlusal stresses during the implant osseointegration timeline, and to stimulate the tissues, a healing crown milled from a resilient hybrid biomaterial is set just a few days after implant placement. Patient satisfaction is thus enhanced by delivering a temporary functional tooth, as a space maintainer, which adapts to the antagonist.

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1st APPOINTMENT

A single minimally invasive surgery is performed and an impression is taken. The Z1 implant replaces the lost root of a tooth, naturally respecting cervical tooth-like anatomy.

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2nd APPOINTMENT

A healing crown made of a hybrid biomaterial enables the patient to leave with a temporary tooth during the osseointegration phase. Tissue healing is thus protected and stimulated.

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3rd APPOINTMENT

The prosthetic treatment plan can now be confirmed and the occlusion accurately recorded.

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4th APPOINTMENT

After osseointegration, the final prosthesis, thus guided and balanced, will be produced and set to the patient.

9 Cliniquement prouvé Clinically proven

Investigation clinique in vivo sur 15 ans de la sécurité et performance des implants Z1
15-year in vivo clinical investigation on the safety and performance of Z1 implants

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*Post-Market Surveillance (PMS) data collected since 2012.

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Polyvalence prothétique
Prosthetic versatility

In implantology today, there is a consensus between preserving the health of the biological width and long-term success of implant treatment.

The patented Tissue Level Z1 implant combines a pure titanium body and a zirconia collar. Titanium is the reference material for its flexural strength and its osseointegration capabilities. Zirconia, on the other hand, offers excellent biocompatibility with fibroblasts. This property is particularly important in implantology, as it promotes epithelial and connective adhesion and proliferation, ensuring the stability of the sealed tissues.

In addition, while protecting the bone tissue from bacterial infiltration, this transgingival zirconia collar also provides aesthetic benefits by mimicking the cervical anatomy of a natural tooth.

Therefore, the Z1 implant gives the advantage of being a biomimetic Tissue Level implant, combining materials optimized for the surrounding tissues, with the aim of preserving and protecting the biological width.

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J. Hu, I. Atsuta, Y. Ayukawa, T. Zhou, I. Narimatsu, and K. Koyano, "Effect of titanium or zirconia implant abutments on epithelial attachments after ultrasonic cleaning," *J. Oral Sci.*, vol. 62, no. 3, pp. 331–334, 2020, doi: 10.2334/josnusd.19-0332.

L. Lago, L. da Silva, I. Martinez-Silva, and B. Rilo, "Radiographic Assessment of Crestal Bone Loss in Tissue-Level Implants Restored by Platform Matching Compared with Bone-Level Implants Restored by Platform Switching: A Randomized, Controlled, Split-Mouth Trial with 3-year Follow-up," *Int. J. Oral Maxillofac. Implants*, vol. 34, no. 1, pp. 179–186, Jan. 2019, doi: 10.11607/jomi.6954.

D.-J. Lee, J.-S. Ryu, M. Shimono, K.-W. Lee, J.-M. Lee, and H.-S. Jung, "Differential Healing Patterns of Mucosal Seal on Zirconia and Titanium Implant," *Front. Physiol.*, vol. 10, no. 796, 2019, doi: 10.3389/fphys.2019.00796.

A. Maltagliati, F. Angiero, S. Zaky, S. Blasi, and A. Ottonello, "Reduction of Bacterial Proliferation by Zirconium Collar in Dental Implants," *Annu. Res. Rev. Biol.*, vol. 23, no. 1, pp. 1–8, Jan. 2018, doi: 10.9734/ARRB/2018/38270.

S. Marconcini, E. Giammarinaro, J. A. Correia, A. Maltagliati, F. Salvado, U. Covani, "Clinical Performance of Titanium-Zirconia Tissue-level Implants in Patients with Well- controlled and Poorly-controlled Type 2 Diabetes: A cohort study with chair-side assessment of oxidative stress," *Oral Implantol.*, vol. 15, no. 1, 2022, [Online]. Available: <https://www.oimplantology.org/oimp/article/view/373>.

C. Pierre, G. Bertrand, I.Pavy, C. Rey, O. Benhamou, C. Roques and C. Combes, "Antibacterial Electrodeposited Copper-Doped Calcium Phosphate Coatings for Dental Implants," *J. Funct. Biomater.*, vol. 14, no. 1, p. 20, Dec. 2022, doi: 10.3390/jfb14010020.

L. Rimondini, L. Cerroni, A. Carrassi, and P. Torricelli, "Bacterial Colonization of Zirconia Ceramic Surfaces: An In Vitro and In Vivo Study," *Int. J. Oral Maxillofac. Implants*, vol. 17, no. 6, pp. 793–8, 2002.

P. Schüpbach, "Interfaces Between Tissues and Ceramics," in *Advanced Ceramics for Dentistry*, Elsevier, 2014, pp. 201–217. doi: 10.1016/B978-0-12-394619-5.00010-9.

S. Shahdad, D. Bosshardt, M. Patel, N. Razaghi, A. Patankar, and M. Rocuzzo, "Benchmark performance of anodized vs. sandblasted implant surfaces in an acute dehiscence type defect animal model," *Clin. Oral Implants Res.*, vol. 33, no. 11, pp. 1135–1146, Nov. 2022, doi: 10.1111/clr.13996.

S. Tetè, F. Mastrangelo, A. Bianchi, V. Zizzari, and A. Scarano, "Collagen Fiber Orientation Around Machined Titanium and Zirconia Dental Implant Necks: An Animal Study," *Int. J. Oral Maxillofac. Implants*, vol. 24, no. 1, pp. 52–8, 2009.

T. T. Vianna et al., "Evaluation of peri-implant marginal tissues around tissue-level and bone-level implants in patients with a history of chronic periodontitis," *J. Clin. Periodontol.*, vol. 45, no. 10, pp. 1255–1265, Oct. 2018, doi: 10.1111/jcpe.12999.

C. Pierre, G. Bertrand, C. Rey, O. Benhamou, and C. Combes, "Calcium phosphate coatings elaborated by the soaking process on titanium dental implants: Surface preparation, processing and physical-chemical characterization," *Dent. Mater.*, vol. 35, no. 2, pp. 25–35, Feb. 2019, doi: 10.1016/j.dental.2018.10.005.

S. Suryavanshi and M. Verma, "Comparative evaluation of the soft tissue response and esthetics of the titanium implants with Zirconium collar," Department of Prosthodontics, Crown and Bridge, Maulana Azad Institute of Dental Sciences, New Delhi, Clinical report, 2014.

📍 24 impasse René Couzinet, 31500 Toulouse - FRANCE
☎ +33 (0)5 62 16 71 00
✉ contact@tbr.dental
🌐 www.tbr.dental

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