



# Easy Concentration & Analysis of PAHs in River Water Using MonoTrap TD for Thermal Desorption

MonoTrap RSC18TD (without graphite carbon) is useful for analysis of compounds with high boiling point such as PAHs. Plain silica MonoTrap RSC18TD is optimal for better desorption.

## Procedure



From left, Benzo[*g,h,i*]perylene, Dibenz[*a,h*]anthracene, Indeno[1,2,3-*cd*]pyrene

Above is the comparison data of PAHs (25ppt) in water after being stirred and introduced to a TD between MonoTrap with graphite carbon and without graphite carbon. MonoTrap RSC18- TD (without graphite carbon) shows better desorption and higher sensitivity.

### GC Conditions

<b>System</b>	: GC/MS-Thermal Desorption (T-Dex II)
<b>Column</b>	: InertCap 5MS Sil 0.25 mm I.D. × 30 m, df = 0.25 μm
<b>Col.Temp.</b>	: 50 °C (3 min) - 5 °C/min - 310 °C
<b>Carrier Gas</b>	: He, 1mL/min (constant flow)
<b>Desorb Temp.</b>	: 250 °C
<b>Time</b>	: 10 min
<b>Flow</b>	: 1 mL/min
<b>Split</b>	: 10 mL/min
<b>Detection</b>	: MS FASST mode

Concentration in pure water


※0.025 – 0.067 ppb

※※0.25 ppb, n5



### 1. Composition of Calibration Curve and Recovery Rate

Compound	Linear range [ng/mL]	r <sup>2</sup>	Recovery(*)	RSD(**)
Acenaphthylene	0.05-5	0.9991	72 %	9 %
Acenaphthene	0.067-6.7	0.9975	81 %	8 %
Fluorene	0.015-1.5	0.9998	90 %	7 %
Phenanthrene	0.011-1.1	0.9996	105 %	5 %
Anthracene	0.003-0.3	0.9999	100 %	4 %
Fluoranthene	0.025-2.5	0.9983	113 %	4 %
Pyrene	0.027-2.7	0.9977	109 %	4 %
Benz[a]anthracene	0.013-1.3	0.9989	113 %	4 %
Chrysene	0.012-1.2	0.9984	108 %	6 %
Benzo[b]fluoranthene	0.013-1.3	0.9990	117 %	5 %
Benzo[k]fluoranthene	0.015-1.5	0.9986	110 %	8 %
Benzo[a]pyrene	0.016-11.6	0.9993	112 %	5 %
Benzo[g,h,i]perylene	0.014-1.4	0.9990	104 %	5 %
Dibenz[a,h]anthracene	0.011-1.1	0.9872	102 %	6 %
Indenno[1,2,3-cd]pyrane	0.012-1.2	0.9954	87 %	4 %

### 2. PAHs Recovered by MonoTrap (Concentration in Water 0.25 ppb)



### 3. PAHs Quantification in River Water



River water was trapped under the fine weather as the sample within a same prefecture and stirred with MonoTrap RSC18-TD. PAHs between 0.3 ppt - 4 ppt were detected.

Selected ion monitoring chromatogram of *Fluoranthene*

GL Sciences disclaims any and all responsibility for any injury or damage which may be caused by this data directly or indirectly. We reserve the right to amend this information or data at any time and without any prior announcement.

#### **GL Sciences, Inc. Japan**

22-1 Nishishinjuku 6-Chome  
Shinjuku-ku, Tokyo,  
163-1130, Japan  
Phone: +81-3-5323-6620  
Fax: +81-3-5323-6621  
Email: [world@glsc.co.jp](mailto:world@glsc.co.jp)  
Web: [www.glsciences.com](http://www.glsciences.com)

#### **GL Sciences B.V.**

De Sleutel 9  
5652 AS Eindhoven  
The Netherlands  
Phone: +31 (0)40 254 95 31  
Email: [info@glsciences.eu](mailto:info@glsciences.eu)  
Web: [www.glsciences.eu](http://www.glsciences.eu)

#### **GL Sciences (ShangHai) Ltd.**

Tower B, Room 2003,  
Far East International Plaza,  
NO,317 Xianxia Road,  
Changning District.  
Shanghai, China P.C. 200032  
Phone: +86 (0)21-6278-2272  
Email: [contact@glsciences.com.cn](mailto:contact@glsciences.com.cn)  
Web: [www.glsciences.com.cn](http://www.glsciences.com.cn)

#### **GL Sciences, Inc. USA**

4733 Torrance Blvd. Suite 255  
Torrance, CA 90503  
Phone: 310-265-4424  
Fax: 310-265-4425  
Email: [info@glsciencesinc.com](mailto:info@glsciencesinc.com)  
Web: [www.glsciencesinc.com](http://www.glsciencesinc.com)

#### **International Distributors**

Visit our Website at:  
<https://www.glsciences.com/company/distributor.html>

