


Detection and determination of 11-nor-9-carboxy- Δ^9 -carboxylic acid (THCA), which is a metabolite of principal psychoactive constituent of cannabis plant, is necessary to confirm cannabis use from urine. Sample pretreatment method suitable for handling low volume and low concentration sample has been desired for the determination because THCA

concentration in urine is in ppb level.

In this note, MonoSpin, which provides high recovery even from sample of less than 100 μL , was used for the sample pretreatment, and the purified solution was injected into LC/MS/MS system. The results showed good linearity, recovery, and reproducibility.

(Y. Yui and S. Ota)


Evaluation of purification by MonoSpin C18-CX (analyzed by HPLC-UV system)



Conditions

Column : Inertsil ODS-3 (5 μm , 150 \times 4.6 mm I.D.)
Eluent: : A) 0.1 % HCOOH in H₂O
 : B) 0.1 % HCOOH in CH₃CN
 A/B = 60/40 - 10 min - 0/100, v/v
Flow rate: Col. : 1.0 mL/min
Temp.: : 40 °C
Detection: : UV 214 nm
Injection Vol.: : 10 μL
Sample: : (\pm)-Carboxy-THC 10 $\mu\text{g}/\text{mL}$

Chemical Structure




(\pm)-11-nor-9-Carboxy- Δ^9 -THC
(THC-COOH)

Structures are created using Chemistry 4-D Draw which is provided by ChemInnovation Software, Inc.

What is MonoSpin ?

MonoSpin is a series of spin columns for solid phase extraction (SPE). Owing to the high permeability of monolithic silica disk packed into the spin column, the procedures, such as conditioning, sample loading, washing, and elution can be carried out only by centrifuging the column. It is also the advantage that the elution volume is only 200 μL .

An Example of Procedures



Elution



Sample Pretreatment using MonoSpin C18-CX

MonoSpin C18-CX has octadecyl group and cation-exchange group on the surface of its silica monolithic support. In this case, acidic buffer was used for sample loading. THC-COOH, which hardly dissociates under acidic condition, could be retained with hydrophobic interaction.

Sample solution	
0.1 % formic acid in water	490 µL
Sample matrix	490 µL
Standard solution *	10 µL
+) Solution of internal standard*	10 µL
	1000 µL

* Their concentration is described in the results.

1. Conditioning (Equilibrating)

Wash with 300 µL of methanol and 300 µL of 0.1 % formic acid in water

Centrifuge**
2 min each

2. Sample loading

Centrifuge

Load 200 µL of sample solution

2 min

3. Wash

Wash with 300 µL of 0.1 % formic acid in water

Centrifuge

2 min

4. Elution

Centrifuge


Elute with 100 µL of 5% ammonia in methanol

1 min

5. HPLC Analysis

** All the centrifugation was operated at 3,000 xg.

Evaluation of Recovery by Standard Addition



Conditions;

System	: LC800 (GL Sciences) API3000 (AB SCIEX)
Column	: InertSustain C18 (2 µm, 100 × 2.1 mm I.D.)
Eluent	: A) 0.1% HCOOH in H ₂ O B) 0.1% HCOOH in CH ₃ CN A/B = 35/65 - 6 min - 5/95 - 4 min - 5/95 - 0.1 min - 35/65 - 5 min - 35/65, v/v
Flow Rate	: 200 µL/min
Col. Temp.	: 40 °C
Detection	: MS/MS (ESI, Nega, MRM)
Injection Vol.	: 5 µL

*** Deuterium labeled compound ((±)-11-nor-9-Carboxy-Δ⁹-THC-D₉, THC-COOH-d₉) was used as internal standard.



Scan Parameters

Analyte	Precursor Ion (m/z)	Product Ion (m/z)
(±)-Carboxy-THC	343.0	299.0
(±)-Carboxy-THC-d ₉	352.1	308.1


Calibration Curves

Calibration curve was plotted 4 times by conducting same experiment on different days. All the calibration curve showed good linear response ($r^2 = 0.9998 \sim 1.000$).


1st day



3rd day



4th day



Recovery, Reproducibility, and Accuracy

From the chromatograms, recovery, reproducibility, and accuracy were calculated for each day. Almost all results were satisfactory.


Recovery (%)	1 st day	2 nd day	3 rd day	4 th day	Inter-day average
0.5 ng/mL	86.8	94.4	87.0	88.5	89.2
1.0 ng/mL	99.0	100.8	101.4	101.7	100.7
10 ng/mL	88.1	91.1	91.6	92.8	90.9

Reproducibility (%)	1 st day	2 nd day	3 rd day	4 th day	Inter-day average
0.5 ng/mL	2.67	0.74	2.17	3.20	2.2
1.0 ng/mL	1.00	2.85	2.79	0.57	1.8
10 ng/mL	1.14	1.05	1.79	0.80	1.2

Accuracy (%)	1 st day	2 nd day	3 rd day	4 th day	Inter-day average
0.5 ng/mL	101.0	105.3	101.0	101.7	102.2
1.0 ng/mL	99.2	97.1	99.5	98.9	98.7
10 ng/mL	99.9	99.9	99.9	100.1	100.1


The series of MonoSpin

MonoSpin C18




Octadecyl group is chemically bonded, and non-polar compounds can be retained because of its hydrophobic interaction. It can be used for extraction or desalting.

MonoSpin SAX




Trimethylaminopropyl group is bonded. It offers strong anion-exchange and weak hydrophobic interaction. It is suitable for extraction of acidic drugs.

MonoSpin NH₂




Aminopropyl group is bonded. It is suitable for extraction of hydrophilic compounds, such as sugar chain.

MonoSpin PBA




Phenylboronic acid is chemically bonded. Compounds containing *cis*-diol group can be retained with high selectivity.

MonoSpin SCX




Propylbenzenesulfonic acid is bonded. It offers strong cation-exchange and hydrophobic interaction. It is suitable for extraction of basic drugs.

MonoSpin TiO




Monolithic silica is coated with titanium dioxide. It is suitable for extraction of phosphate-containing compounds.

MonoSpin Amide




Amide group is bonded. Wide variety of hydrophilic compounds, such as sugar chain, can be retained and extracted in HILIC mode.

MonoSpin CBA



Carboxyl group is bonded to silica monolith. It is suitable for extraction of basic compounds through weak cation-exchange interaction.


MonoSpin C18-CX



Functional groups of both C18 and SCX are bonded. C18-CX is often superior to C18 or SCX in purification and clean-up of basic drugs in serum and urine.

(used in this note)

MonoSpin Trypsin



Trypsin is immobilized on the surface of silica monolithic support. It enables rapid protein-digesting.

GL Sciences disclaims any and all responsibility for any injury or damage which may be caused by this data directly or indirectly. We reserve the right to amend this information or data at any time and without any prior announcement.

GL Sciences, Inc. Japan

22-1 Nishishinjuku 6-Chome
Shinjuku-ku, Tokyo,
163-1130, Japan
Phone: +81-3-5323-6620
Fax: +81-3-5323-6621
Email: world@glsc.co.jp
Web: www.glsciences.com

GL Sciences B.V.

De Sleutel 9
5652 AS Eindhoven
The Netherlands
Phone: +31 (0)40 254 95 31
Email: info@glsciences.eu
Web: www.glsciences.eu

GL Sciences (ShangHai) Ltd.

Tower B, Room 2003,
Far East International Plaza,
NO,317 Xianxia Road,
Changning District.
Shanghai, China P.C. 200032
Phone: +86 (0)21-6278-2272
Email: contact@glsciences.com.cn
Web: www.glsciences.com.cn

GL Sciences, Inc. USA

4733 Torrance Blvd. Suite 255
Torrance, CA 90503
Phone: 310-265-4424
Fax: 310-265-4425
Email: info@glsciencesinc.com
Web: www.glsciencesinc.com

International Distributors

Visit our Website at:

<https://www.glsciences.com/company/distributor.html>