


## Analysis of Amines



### Conditions

<b>System</b>	: GC-4000 Plus system (GL Sciences Inc.)	<b>Analyte :</b>	
<b>Column</b>	: InertCap for Amines (GL Sciences Inc.) 0.32 mm I.D. x 30 m	1. Monomethylamine	500 ppm (v/v)
<b>Column Cat. No.</b>	: 1010-69249	2. Dimethylamine	500 ppm (v/v)
<b>Column Temp.</b>	: 35 °C(3 min) - 15 °C/min - 260 °C(22 min)	3. Trimethylamine	500 ppm (v/v)
<b>Carrier Gas</b>	: He 100 kPa	4. Triethylamine	0.5 mg/mL
<b>Injection</b>	: Split 60 mL/min 260 °C	5. Tri- <i>n</i> -butylamine	0.5 mg/mL
<b>Detection</b>	: FID Auto Range 260 °C	6. Dodecylamine	0.5 mg/mL
<b>Injection Vol.</b>	: 1 µL	7. Hexadecylamine	
<b>Sample</b>	: Standard	8. Octadecylamine	
		Sol. in 2-Propanol	