

How to Use Preparative HPLC - Part 2


Scaling up from Analytical HPLC

It is not easy to find out optimal condition for preparative HPLC. Not only large volume of solvent but also substantial amount of precious sample may be required for the evaluation of separation conditions, particularly in preparative HPLC. Consequently, we recommend that the evaluation should be carried out using analytical column (4.6 mm I.D.) in the beginning. Condition for preparative HPLC can be investigated efficiently by using analytical column packed with the same gel as in preparative HPLC column.

In this note, Inertsil ODS-3 was taken as an example, and how to scaling up from analytical column to preparative column is described.

(K. Kanno)

Conditions of preparative HPLC are generally optimized as the following flowchart



<What is important to scale up>

It is important to calculate ratio of cross-sectional area of preparative column to that of analytical column. The ratio can be used as follows;

1) Increase flow rate in proportion to cross-sectional area of column

2) Increase loading amount (injection volume) in proportion to cross-sectional area of column


Inner diameter	4.6 mm	10 mm	14 mm	20 mm	30 mm
Cross-sectional area	16.6 mm ²	78.5 mm ²	154 mm ²	314 mm ²	707 mm ²
Ratio of the area※ ¹	1	4.7	9.3	18.9	42.5

※¹ Ratio of cross-sectional to that of 4.6-mm I.D. column was calculated.

1) Increase flow rate in proportion to cross-sectional area of column

The figure shown right indicates relation between linear velocity of mobile phase and height equivalent to theoretical plate (HETP) obtained with three columns packed with 10 μm particles.

Optimum flow rate, at which the lowest HETP is obtained, is 3.0 cm/min (0.5 mm/sec) for all the columns. Therefore, it can be said that flow rate should be changed to maintain optimum linear velocity of 3.0 cm/min in case of scaling up from 10 μm particle packed analytical column to 10 μm particle packed preparative one. It is important that particle size of the two columns is same because optimum flow rate changes also depending on particle size.




2) Increase loading amount (injection volume) in proportion to cross-sectional area of column

The figure shown below indicates relation between loading amount and number of theoretical plates (N). Three columns with different inner diameters were used and compared. For example, maximum loading amount for each column at which N above 2000 can be obtained is follows;


6 mm I.D. approx. 1 mg
20 mm I.D. approx. 10 mg
50 mm I.D. approx. 70 mg

Since the maximum loading amount is proportional to cross-sectional area of column, it can be said that similar separation should be achieved with preparative column as with analytical column by increasing injection volume in proportion to cross-sectional area.



<Parameters to be changed for scaling up>

In case of scaling up from a 4.6 mm I.D. analytical column to a 20 mm I.D. preparative column, cross-sectional area of the column is approximately 19 times enlarged. Therefore, scaling up can be achieved by increasing flow rate and loading amount (injection volume) 19 times. **Red letters** represent parameters to be changed for scaling up.

Column	: Inertsil ODS-3 (10 μ m, 4.6 mm I.D. \times 250 mm)	Scale up	Column	: Inertsil ODS-3 (10 μ m, 20 mm I.D. \times 250 mm)
Eluent	: A) CH ₃ CN B) H ₂ O A/B = 40/60, v/v		Eluent	: A) CH ₃ CN B) H ₂ O A/B = 40/60, v/v
Flow rate	: 500 μL/min	Increase flow rate and injection volume 19 times.	Flow rate	: 9.5 mL/min (9500 μ L/min)
Column Temp.	: 40 $^{\circ}$ C	Keep other parameters unchanged	Column Temp.	: 40 $^{\circ}$ C
Detection	: UV 270 nm		Detection	: UV 270 nm
Injection Vol.*	: 500 μL		Injection Vol.*	: 9.5 mL (9500 μL)


* Concentration of the sample solution is same

<An example of scaling up>

Chromatograms before and after scaling up are shown below


HPLC conditions

Column : Inertsil ODS-3
(5 μ m, 250 mm length)
Column I.D. : Discrbed in each chromatogram
Eluent : A) CH₃CN B) 0.1 % TFA
A/B = 40/60, v/v
Flow rate : Discrbed in each chromatogram
Column Temp. : 40 $^{\circ}$ C
Detection : UV 270 nm
Sample : Commercially available licorice extract (2.0 mg/mL)




Column I.D. : 4.6 mm
Flow rate : 1.0 mL/min
Injection Vol. : 160 μ L
Loading amount : 0.32 mg

Loading amount was 18.9 times increased.



Column I.D. : 20 mm
Flow rate : 18.9 mL/min
Injection Vol. : 3.0 mL
Loading amount : 6.0 mg



GL Sciences disclaims any and all responsibility for any injury or damage which may be caused by this data directly or indirectly. We reserve the right to amend this information or data at any time and without any prior announcement.

GL Sciences, Inc. Japan

22-1 Nishishinjuku 6-Chome
Shinjuku-ku, Tokyo,
163-1130, Japan
Phone: +81-3-5323-6620
Fax: +81-3-5323-6621
Email: world@glsc.co.jp
Web: www.glsciences.com

GL Sciences B.V.

De Sleutel 9
5652 AS Eindhoven
The Netherlands
Phone: +31 (0)40 254 95 31
Email: info@glsciences.eu
Web: www.glsciences.eu

GL Sciences (ShangHai) Ltd.

Tower B, Room 2003,
Far East International Plaza,
NO,317 Xianxia Road,
Changning District.
Shanghai, China P.C. 200032
Phone: +86 (0)21-6278-2272
Email: contact@glsciences.com.cn
Web: www.glsciences.com.cn

GL Sciences, Inc. USA

4733 Torrance Blvd. Suite 255
Torrance, CA 90503
Phone: 310-265-4424
Fax: 310-265-4425
Email: info@glsciencesinc.com
Web: www.glsciencesinc.com

International Distributors

Visit our Website at:

<https://www.glsciences.com/company/distributor.html>