


### Analysis of Tranexamic acid

(Under the Condition of the Japanese Pharmacopoeia, Tranexamic acid Capsules)

Data No. LB453-0919



#### Conditions

**System** : GL7700 HPLC system  
**Column** : Inertsil ODS-4  
(5  $\mu$  m, 150 x 4.6 mm I.D.)  
**Column Cat. No.** : 5020-03945  
**Eluent** : A) CH<sub>3</sub>OH  
B) Buffer\*  
A/B = 40/60, v/v  
**Flow Rate** : 0.9 mL/min  
**Col. Temp.** : 25 °C  
**Detection** : UV 220 nm (PD7752 PDA Detector)  
**Injection Vol.** : 10  $\mu$  L  
**Sample** : Standard

\*Dissolve 11.0 g of anhydrous sodium dihydrogen phosphate in 500 mL of water, and add 10 mL of triethylamine and 1.4 g of sodium lauryl sulfate. Adjust pH 2.5 with phosphoric acid, add water to make 600 mL.

#### Analyte:

1. Tranexamic acid 280 mg/L

Theoretical plate number : 6,899 ( $\geq$  4,000)  
Tailing factor : 1.07 ( $\leq$  2.0)  
RSD of the peak area (n=6) : 0.25 ( $\leq$  2.0)

#### 【NOTE】

- 1) Fully equilibrate the column prior to the analysis. Fully equilibrate the column with eluent for at least 24 hrs at 1 mL/min.
- 2) Prepare the eluent at time of use, otherwise the retention time may shift.