


Analysis of Haloperidol

(Under the Condition of the Japanese Pharmacopoeia 17th edition)

Data No. LB431-0894


System performance

1-a. Haloperidol 0.5 mg/L




Standard solution

1-b. Haloperidol 2.5 mg/L



Sample solution

1-c. Haloperidol 500 mg/L



Conditions

System : GL7700 HPLC system
Column : Inertsil ODS-4
 (5 μm, 150 x 4.6 mm I.D.)
Column Cat. No. : 5020-03945
Eluent : Buffer*
Flow Rate : 1.0 mL/min
Col. Temp. : 40 °C
Detection : UV 220 nm (UV7750 UV Detector)
Injection Vol. : 10 μL
Sample : Standard

Analyte

1. Haloperidol

The peak area ratio of 1-a to 1-b : (15 ≤) 20.0 (≤ 25)

Theoretical plates of 1-b : 9,176 (≥ 4,000)

Tailing factor of 1-b : 1.06 (≤ 2.0)

RSD of the
 1-b area (%) (n=6) : 0.55 (≤ 2.0)

*Dissolve 2.95 g of trisodium citrate dihydrate in 900 mL of water.

Adjust to pH 3.5 with diluted hydrochloric acid.

Add water to make 1,000 mL. (Solution A)

Add 700 mL of CH₃OH to 300 mL of solution A. (Solution B)

Add 1.0 g of sodium dodecyl sulfate to Solution B.