

Bioseparation Solution

Bio LC Column

ProteoSil / MonoSelect



Peptide Mapping
Nucleic acids / oligonucleotides
Monoclonal Antibodies



Introduction

GL Sciences Bio LC Columns are HPLC solutions specifically designed for bio molecules proteomics such as protein, peptide and nucleic acid species.

All columns are available with bio-inert PEEK and stainless steel hardware.


Packing materials are made of high-purity silica with pore sizes of 100Å, 200Å, and 300Å, and are available in reversed phase, HILIC, and size exclusion columns. Low lot-to-lot variation and consistent quality also make it ideal for LC/MS applications.

| | Target | MW | Column | Phase | Particle size (um) | Pore size (Å) |
|----------------|--|---------------|-------------------|------------------------|--------------------|---------------|
| Reversed Phase | Proteins / Peptides Oligonucleotides / Nucleic Acids | >20,000 | ProteoSil 300-C18 | C18 | 5 | 300 |
| | | | ProteoSil 300-C8 | C8 | 5 | 300 |
| | | 5,000-20,000 | ProteoSil 200-C18 | C18 | 1.9, 3, 5 | 200 |
| | | <5,000 | ProteoSil 100-C18 | C18 | 1.9, 3, 5 | 100 |
| | Proteins / Peptides Monoclonal Antibodies Oligonucleotides / Nucleic Acids | >20,000 | ProteoSil 300-C4 | C4 | 5 | 300 |
| | Proteins, Monoclonal Antibodies Antibody-Drug Conjugate(ADC) Sub unit | >100,000 | MonoSelect RP-mAb | Phenyl | Monolith | 600 |
| HILIC | Proteins / Peptides Monoclonal Antibodies | — | ProteoSil HILIC | Amide | 1.9, 3, 5 | 100 |
| SEC | Proteins / Peptides Monoclonal Antibodies Oligonucleotides / Nucleic Acids | 5,000-600,000 | ProteoSil 300-SEC | DIOL | 5 | 300 |
| | | <5,000 | ProteoSil 100-SEC | DIOL | 5 | 100 |
| SEC + RP | LNP, Exosome | — | MonoSelect nPEC | Hydrophilic Polymer | Monolith | 110 |


ProteoSil Pore Size Range

ProteoSil HPLC columns feature precisely controlled pore sizes in their packing material, making them ideal for the analysis of biopharmaceuticals. Specifically, the 200Å pore size is optimized for the analysis of compounds with molecular weights ranging from a few kDa to several tens of kDa, delivering optimal performance for the analysis of peptides and oligonucleotides.

Packing material pore size and analyte molecular weight range



Comparison of 200Å and 175Å columns for oligonucleotide analysis



Column Hardware

The column hardware can be selected from stainless steel and Bio-Inert PEEK. Innovative PEEK-lined stainless steel has increased the maximum operating pressure.



Steel-Coated-PEEK
(UHPLC PEEK)




Bioinert PEEK



Stainless Steel

Guard Column

The guard column is installed between injector and analysis column, it is mainly used to protect the analysis column. There are 3 types guard column



Cartridge Guard Column E



Cartridge Guard Column Ei
(Non-metal Type)





Guard Columns for UHPLC


Quality

To ensure that our HPLC columns consistently maintain the same high quality and are reliably supplied worldwide, our production facility holds ISO 9001 and ISO 14001 certifications. We carry out all processes, including matrix synthesis, chemical treatment, column packing, and column quality inspections in-house. We continuously evolve based on our accumulated expertise, enabling us to offer superior HPLC columns.

In addition to the chemical modification of functional groups and end-cap processing, we also synthesize silica gel, a critical component for column performance.



Reliable Lot-to-Lot Consistency



Product Line

ProteoSil 300-C18

Designed for advanced protein analysis, the ProteoSil 300-C18 column offers the ideal capabilities when enhanced retention and sample load are critical for C4 or C8 columns.

ProteoSil 300-C8

Optimized for proteins and peptides, the ProteoSil 300-C8 ensures unique selectivity and hydrophobicity compared to shorter-chain C4 phases, suitable for protein and hydrophobic peptide separations.

ProteoSil 200-C18

Ideal for mid-sized molecules, the ProteoSil 200-C18 column allows for separation across 100Å pore columns and is easy to clean for faster analysis.

ProteoSil 200-C8

Specialized for the analysis of mid-sized molecules, the ProteoSil 200-C8 shows superior retention for analytes that are strongly retained by C18 columns.

ProteoSil 100-C18

Offering high inertness and durability, the ProteoSil 100-C18 is perfect for analyzing low-molecular-weight compounds, enhancing sensitivity in LC-MS applications.

ProteoSil 100-C8

Suitable for low-molecular-weight peptides, the ProteoSil 100-C8 provides high inertness and durability with lower retention compared to C18 columns.

ProteoSil 300-C4

Tailored for the analysis of proteins and peptides, especially hydrophobic peptides, the ProteoSil 300-C4 facilitates shorter retention times and is ideal for proteins with moderate retention needs.

ProteoSil HILIC

Modified with amide groups, the ProteoSil HILIC column excels in the separation of highly hydrophilic compounds, peptides, glycans, and oligonucleotides.

ProteoSil 300-SEC

Equipped with dihydroxypropyl groups bonded to silica gel, the ProteoSil 300-SEC column is suited for the analysis of large molecules and can be used as a dial column.

ProteoSil 100-SEC

Constructed for aqueous size-exclusion chromatography, the ProteoSil 100-SEC can handle higher operating pressures and is suitable for series column configurations.

MonoSelect RP-mAb

Featuring a silica monolith structure with phenyl groups, the MonoSelect RP-mAb column is optimized for monoclonal antibodies analysis and high-sensitivity LC/MS applications.

MonoSelect nPEC

Designed for rapid separation of nanoparticles and free low-molecular-weight compounds, the MonoSelect nPEC column is crucial for assessing drug encapsulation rates in liposome formulations.




Reversed Phase

C18 Columns are a recommended for peptide and protein analysis.


It uses a packing material with a pore size of 200Å, making it ideal in the determination of low to high molecular weight compounds (up to several 10kDa). The use of extremely non-specific binding to packing materials and metal-free column hardware enables analyte sharp peaks even for adsorbable analytes.

Analysis of Peptides



Conditions

| | | |
|----------------|--|--------------------------|
| Column | : ProteoSil 200-C18 UHPLC PEEK (1.9 μm, 100 × 2.1 mm I.D) | Sample : |
| Eluent | : A) 0.1% HCOOH in H ₂ O B) 0.1% HCOOH in CH ₃ CN | 1. Gly-Tyr |
| | A/B = 95/5 - 0.5 min - 70/30 - 2.5 min - 60/40 | 2. Val-Tyr-Val |
| | - 0.5 min - 60/40 - 0.01/min - 95/5 - 6.49 min - 95/5, v/v | 3. Angiotensin II |
| Flow Rate | : 0.3 mL/min | 4. Methionine enkephalin |
| Col. Temp. | : 40°C | 5. Leucine enkephalin |
| Detection | : UV 280 nm | (50 mg/mL each) |
| Injection Vol. | : 5 μL | |




Conditions

| | |
|----------------|--|
| Column | : ProteoSil 200-C18 UHPLC PEEK (1.9 μm, 100 × 2.1 mm I.D) |
| Eluent | : A) 0.1 % Triethylamine in H ₂ O (pH 6.3, CH ₃ COOH) B) Eluent A/CH ₃ CN = 50/50, v/v |
| | A/B = 83/17 - 4 min - 80/20 - 0.1 min - 83/17 - 5.9 min - 83/17, v/v |
| Flow Rate | : 0.4 mL/min |
| Col. Temp. | : 40°C |
| Detection | : UV 260 nm |
| Injection Vol. | : 10 μL |
| Sample | : 1. 5' - GTT ACA GAA TCT GAC AAG CCT AAT ACG - 3' (27 mer) 2. 5' - GTT ACA GAA TCT GCC AAG CCT AAT ACG - 3' (27 mer) 3. 5' - GTT ACA GAA TCT GTC AAG CCT AAT ACG - 3' (27 mer) (300 pmol/L each) |

Analysis of Monoclonal Antibodies (mAbs)


MonoSlect RP-mAb is a HPLC column Specialized for Monoclonal Antibody Analysis. Monolithic silica consists of co-continuous through- pores and skeletons which have mesopores. The large surface area and high permeability of this structure enables “strong retentivity and low pressure”. Broad peaks are obtained with conventional HPLC columns when the analytes are huge proteins larger than 10 nm such as antibodies. “Sharp peaks of mAb” can be obtained with MonoSelect RP-mAb because the mesopores are designed to be 60 nm, which is suitable for mAb analysis. The analytical time is so rapidly.



Reduction process

Conditions

| | |
|----------------|--|
| Column | : MonoSelect RP-mAb |
| Eluent | : A) 0.075 % HCOOH + 0.025 % TFA in H ₂ O B) 0.075 % HCOOH + 0.025 % TFA in CH ₃ CN |
| | A/B = 80/20-(5min)-40/60,v/v |
| Flow Rate | : 0.3 mL/min |
| Column Temp. | : 80 °C |
| Detection | : 214 nm |
| Injection Vol. | : 10 μL |



Papain digestion process


Conditions

| | |
|----------------|--|
| Column | : MonoSelect RP-mAb |
| Eluent | : A) 0.075 % HCOOH + 0.025 % TFA in H ₂ O B) 0.075 % HCOOH + 0.025 % TFA in CH ₃ CN |
| | A/B = 85/15-(5min)-60/40,v/v |
| Flow Rate | : 0.3 mL/min |
| Column Temp. | : 80 °C |
| Detection | : 214 nm |
| Injection Vol. | : 10 μL |

HILIC (Hydrophilic Interaction Liquid Chromatography) Columns

HILIC (Hydrophilic Interaction Liquid Chromatography) mode is used in the analysis of biomolecules. It has a hydrophilic stationary phase and water-based solvent as the mobile phase, making it ideal for separating and analyzing hydrophilic compounds, such as polar biomolecules and glycans. HILIC mode offers high sensitivity analysis when combined with mass spectrometry, making it suitable for detecting trace amounts of biomolecules and specific analyses. It is particularly useful for glycan analysis and can serve as an alternative to ion exchange chromatography in some cases.

Analysis of Peptides



Conditions

Column : ProteoSil HILIC
(1.9 μm , 150 \times 2.1 mm I.D.)
Eluent : A) CH₃CN
B) 10 mM HCOONH₄ + 0.1% HCOOH in H₂O

| Time (min) | A (vol%) | B (vol%) |
|------------|----------|----------|
| 0 | 80 | 20 |
| 15 | 50 | 50 |

Flow Rate : 0.2 mL/min
Col. Temp. : 40 °C
Detection : UV 254 nm
Injection Vol. : 1 μL
Sample : Reference Standard


Analyte :

1. Leucine-Enkephalin 500 mg/L
(Tyr-Gly-Gly-Phe-Leu)
2. Methionine-Enkephalin 500 mg/L
(Tyr-Gly-Gly-Phe-Met)
3. Angiotensin II (Human) 500 mg/L
(Asp-Arg-Val-Tyr-Ile-His-Pro-Phe)
4. Val-Tyr-Val 500 mg/L
5. Gly-Tyr 500 mg/L


Size Exclusion Chromatography (SEC) Columns

The ProteoSil SEC column is equipped with a dihydroxypropyl group bonded to silica gel, making it suitable for the analysis of large biomolecules. In SEC, larger analytes such as proteins elute from the column first, while the smallest molecules that can access the pores elute later from the column. This column is available in variations with pore sizes of 300Å and 100Å.


Calibration Curve and Analysis of Proteins



| MW | | Conditions | |
|-----------------------|---------|----------------|--|
| 1. Thyroglobulin | 660,000 | Column | : ProteoSil 300-SEC (5 μm, 250 × 7.6 mm I.D.) |
| 2. Ferritin | 444,000 | Eluent | : 0.2 M NaCl, 0.1 M Na ₂ HPO ₄ (pH 6.7, Na ₂ HPO ₄) |
| 3. γ-Globulins | 160,000 | Flow Rate | : 1.0 mL/min |
| 4. Phosphorylase | 97,400 | Col. Temp. | : 35 °C |
| 5. Transferrin, human | 80,000 | Detection | : UV 280 nm |
| 6. BSA | 66,000 | Injection Vol. | : 10 μL |
| 7. Ovalbumin | 45,000 | | |
| 8. Trypsinogen | 25,000 | | |
| 9. Lysozyme | 14,300 | | |
| 10. Insulin | 6,000 | | |
| 11. Insulin chain B | 3,500 | | |



Analysis of doxorubicin




| Conditions | |
|----------------|--|
| Column | : MonoSelect nPEC (50 × 3.0 mm I.D.) |
| Eluent | : A) 50 mM Na ₂ SO ₄ in 20 mM Acetate Buffer (pH4.6) B) CH ₃ OH A/B = 95/5, v/v |
| Flow Rate | : 1.0 mL/min |
| Column Temp. | : 30 °C |
| Detection | : 254 nm |
| Injection Vol. | : 5 μL |

*This product is developed in collaboration with Dr. Kato, School of Pharmacy, Showa University, and Eisai Co., Ltd.


Comparison Data -MonoSelect mAb-

In protein analysis, especially in low-concentration analysis, the peak shape may deteriorate. In this case, the pore size of sorbent is the most important. The MonoSelect RP-mAb is used monolithic gel instead of particle silica gel. The mesopore of monolithic gel are designed to be optimal for retaining mAb at 60nm and are also treated with an optimal inert treatment to prevent non-specific adsorption on the surface. This design is to provide sharp peak shape even in low-concentration analysis.


MonoSelect RP-mAb Pore Size: 60 nm
2.1 × 20 mm



Brand W Pore Size: 45 nm
2.7 μm, 2.1 × 50 mm



Brand A Pore Size: 100 nm
5 μm, 2.1 × 50 mm



Initial Pressure: **1.7MPa**

4.5 MPa

10.5 MPa

| | |
|----------------|---|
| Eluate | : A) 0.075 % HCOOH + 0.025 % CF ₃ COOH in H ₂ O B) 0.075 % HCOOH + 0.025 % CF ₃ COOH in CH ₃ CN A/B = 95/5-(5min)-10/90 |
| Injection Vol. | : 1μL |
| Oven Temp. | : 80 °C |
| Analyte | : 1.4 mg/mL mAb |

Ordering Information -Reversed Phase-

Maximum Operating Pressure

| ID (mm) | Particle Size (µm) | Maximum Operating Pressure |
|-----------|--------------------|----------------------------|
| 2.1 - 3.0 | 1.9 | 80 MPa (800 bar) |
| 2.1 - 4.6 | 3 HP | 50 MPa (500 bar) |
| 1.0 - 4.6 | 3, 5 | 20 MPa (200 bar) |

ProteoSil 100-C18

| Hardware | Particle Size (µm) | Length (mm) | ID (mm) | |
|----------|-------------------------|-------------|------------|------------|
| | | | 2.1 | 4.6 |
| SS | 1.9 | 50 | 5020-42288 | - |
| | | 75 | 5020-42292 | - |
| | | 100 | 5020-42289 | - |
| | | 150 | 5020-42290 | - |
| | 3 50 MPa (500bar) | 50 | 5020-42278 | 5020-42283 |
| | | 75 | 5020-42279 | 5020-42284 |
| | | 100 | 5020-42280 | 5020-42285 |
| | | 150 | 5020-42281 | 5020-42286 |
| | | 250 | 5020-42282 | 5020-42287 |
| | 5 | 50 | 5020-42266 | 5020-42271 |
| | | 75 | 5020-42267 | 5020-42272 |
| | | 100 | 5020-42268 | 5020-42273 |
| | | 150 | 5020-42269 | 5020-42274 |
| | | 250 | 5020-42270 | 5020-42275 |
| | UHPLC PEEK | 1.9 | 50 | 5020-42264 |
| 100 | | | 5020-42265 | 5020-42262 |
| 150 | | | - | 5020-42263 |
| 3 | | 50 | 5020-42257 | 5020-42253 |
| | | 100 | 5020-42258 | 5020-42254 |
| 150 | 5020-42259 | 5020-42255 | | |
| 250 | 5020-42260 | 5020-42256 | | |
| PEEK | 5 | 50 | 5020-42249 | 5020-42245 |
| | | 100 | 5020-42250 | 5020-42246 |
| | | 150 | 5020-42251 | 5020-42247 |
| | | 250 | 5020-42252 | 5020-42248 |

ProteoSil 200-C18

| Hardware | Particle Size (µm) | Length (mm) | ID (mm) | |
|----------|-------------------------|-------------|------------|------------|
| | | | 2.1 | 4.6 |
| SS | 1.9 | 50 | 5020-42214 | - |
| | | 75 | 5020-42215 | - |
| | | 100 | 5020-42216 | - |
| | | 150 | 5020-42217 | - |
| | 3 50 MPa (500bar) | 50 | 5020-42204 | 5020-42209 |
| | | 75 | 5020-42205 | 5020-42210 |
| | | 100 | 5020-42206 | 5020-42211 |
| | | 150 | 5020-42207 | 5020-42212 |
| | | 250 | 5020-42208 | 5020-42213 |
| | 5 | 50 | 5020-42191 | 5020-42196 |
| | | 75 | 5020-42192 | 5020-42197 |
| | | 100 | 5020-42193 | 5020-42198 |
| | | 150 | 5020-42194 | 5020-42199 |
| | | 250 | 5020-42195 | 5020-42200 |
| | UHPLC PEEK | 1.9 | 50 | 5020-42178 |
| 100 | | | 5020-42179 | 5020-42181 |
| 150 | | | - | 5020-42182 |
| PEEK | 5 | 50 | 5020-42183 | 5020-42187 |
| | | 100 | 5020-42184 | 5020-42188 |
| | | 150 | 5020-42185 | 5020-42189 |
| | | 250 | 5020-42186 | 5020-42190 |

ProteoSil 300-C18

| Hardware | Particle Size (µm) | Length (mm) | ID (mm) | |
|----------|--------------------|-------------|------------|------------|
| | | | 2.1 | 4.6 |
| SS | 5 | 50 | 5020-42110 | 5020-42115 |
| | | 75 | 5020-42111 | 5020-42116 |
| | | 100 | 5020-42112 | 5020-42117 |
| | | 150 | 5020-42113 | 5020-42118 |
| | | 250 | 5020-42114 | 5020-42119 |

Ordering Information -Reversed Phase-

ProteoSil 200-C8

| Hardware | Particle Size (µm) | Length (mm) | ID (mm) | |
|------------|-------------------------|-------------|------------|------------|
| | | | 2.1 | 4.6 |
| SS | 3 50 MPa (500bar) | 75 | 5020-42301 | 5020-42302 |
| | | 50 | 5020-42303 | 5020-42307 |
| | | 75 | 5020-42312 | 5020-42313 |
| | | 100 | 5020-42304 | 5020-42308 |
| | | 150 | 5020-42305 | 5020-42309 |
| | | 250 | 5020-42306 | 5020-42310 |
| UHPLC PEEK | 3 | 50 | 5020-42293 | 5020-42297 |
| | | 100 | 5020-42294 | 5020-42298 |
| | | 150 | 5020-42295 | 5020-42299 |
| | | 250 | 5020-42296 | 5020-42300 |

ProteoSil 300-C8

| Hardware | Particle Size (µm) | Length (mm) | ID (mm) | |
|----------|--------------------|-------------|------------|------------|
| | | | 2.1 | 4.6 |
| SS | 5 | 50 | 5020-42100 | 5020-42105 |
| | | 75 | 5020-42101 | 5020-42106 |
| | | 100 | 5020-42102 | 5020-42107 |
| | | 150 | 5020-42103 | 5020-42108 |
| | | 250 | 5020-42104 | 5020-42109 |

ProteoSil 300-C4

| Hardware | Particle Size (µm) | Length (mm) | ID (mm) | |
|----------|--------------------|-------------|------------|------------|
| | | | 2.1 | 4.6 |
| SS | 5 | 50 | 5020-42120 | 5020-42125 |
| | | 75 | 5020-42121 | 5020-42126 |
| | | 100 | 5020-42122 | 5020-42127 |
| | | 150 | 5020-42123 | 5020-42128 |
| | | 250 | 5020-42124 | 5020-42129 |

MonoSelect RP-mAb set (Holder+Cartridge)

| Item | ID (mm) | Length (mm) | Cat.No. |
|------------------------------|---------|-------------|------------|
| MonoSelect RP-mAb Holder Set | 2.1 | 20 | 5020-10818 |

MonoSelect RP-mAb Cartridge

| Item | ID (mm) | Length (mm) | Cat.No. |
|-----------------------------|---------|-------------|------------|
| MonoSelect RP-mAb Cartridge | 2.1 | 20 | 5020-10819 |



MonoSelect RP-mAb

Cartridge Holder

| Item | Length of the Cartridge Applicable | Cat.No. |
|-----------------------------|------------------------------------|------------|
| MonoSelect Cartridge Holder | 20 mm | 5020-10815 |

Ordering Information -Size Extraction Chromatography (SEC)-

ProteoSil 100-SEC

| Hardware | Particle Size (μm) | Length (mm) | ID (mm) | | |
|----------|--------------------|-------------|------------|------------|------------|
| | | | 4.6 | 6 | 7.6 |
| SS | 5 | 50 | 5020-42315 | 5020-42320 | 5020-42325 |
| | | 75 | 5020-42316 | - | - |
| | | 100 | 5020-42317 | 5020-42321 | 5020-42326 |
| | | 150 | 5020-42318 | 5020-42322 | 5020-42327 |
| | | 250 | 5020-42319 | 5020-42323 | 5020-42328 |

ProteoSil 300-SEC

| Hardware | Particle Size (μm) | Length (mm) | ID (mm) | | | |
|----------|--------------------|-------------|------------|------------|------------|------------|
| | | | 2.1 | 4.6 | 6 | 7.6 |
| SS | 5 | 50 | - | 5020-42134 | 5020-42139 | 5020-42144 |
| | | 75 | 5020-42130 | 5020-42135 | - | - |
| | | 100 | 5020-42131 | 5020-42136 | 5020-42140 | 5020-42145 |
| | | 150 | 5020-42132 | 5020-42137 | 5020-42141 | 5020-42146 |
| | | 250 | 5020-42133 | 5020-42138 | 5020-42142 | 5020-42147 |
| PEEK | 5 | 150 | - | 5020-42168 | - | - |
| | | 250 | - | 5020-42169 | - | - |

MonoSelect nPEC set (Holder+Cartridge)


| Item | ID (mm) | Length (mm) | Cat.No. |
|----------------------|---------|-------------|------------|
| MonoSelect nPEC Sets | 3.0 | 50 | 5020-10816 |



DMA-nPEC

MonoSelect nPEC Cartridge

| Item | ID (mm) | Length (mm) | Cat.No. |
|---------------------------|---------|-------------|------------|
| MonoSelect nPEC Cartridge | 3.0 | 50 | 5020-10817 |



DMA-nPEC cartridge

Cartridge Packing

| Item | Qty. | Cat.No. |
|-------------------------|-------|------------|
| MonoSelect nPEC Packing | 6 pcs | 5020-10880 |



DMA-nPEC packing

Ordering Information - HILIC -

ProteoSil HILIC

| Hardware | Particle Size (μm) | Length (mm) | ID (mm) | |
|----------|--------------------|-------------|------------|------------|
| | | | 2.1 | 4.6 |
| SS | 1.9 | 75 | 5020-42238 | - |
| | | 100 | 5020-42239 | - |
| | | 150 | 5020-42240 | - |
| | 3 | 50 | 5020-42233 | - |
| | | 75 | 5020-42234 | - |
| | | 100 | 5020-42235 | - |
| | | 150 | 5020-42236 | - |
| | 5 | 250 | 5020-42237 | - |
| | | 50 | 5020-42220 | 5020-42225 |
| | | 75 | 5020-42221 | 5020-42226 |
| | | 100 | 5020-42222 | 5020-42227 |
| | | 150 | 5020-42223 | 5020-42228 |
| | 250 | 5020-42224 | 5020-42229 | |

| Hardware | Particle Size (μm) | Length (mm) | ID (mm) | |
|------------|--------------------|-------------|------------|------------|
| | | | 2.1 | 4.6 |
| PEEK | 5 | 50 | 5020-42170 | 5020-42174 |
| | | 100 | 5020-42171 | 5020-42175 |
| | | 150 | 5020-42172 | 5020-42176 |
| | | 250 | 5020-42173 | 5020-42177 |
| UHPLC PEEK | 1.9 | 50 | 5020-42242 | - |
| | | 100 | 5020-42243 | - |
| | | 150 | 5020-42244 | - |
| | 3 | 50 | 5020-42160 | 5020-42164 |
| | | 100 | 5020-42161 | 5020-42165 |
| | | 150 | 5020-42162 | 5020-42166 |
| | | 250 | 5020-42163 | 5020-42167 |

Ordering Information -Guard Column-

Replacement Guard Cartridge for ProteoSil UHPLC

| Hardware | Particle Size (µm) | Length (mm) | Phase | ID (mm) | |
|----------|--------------------|-------------|---------|------------|--|
| | | | | 2.1 | |
| SS | 1.9 | 10 | 100-C18 | 5020-42291 | |
| | | | 200-C18 | 5020-42219 | |
| | | | HILIC | 5020-42241 | |

Cartridge Holder for UHPLC

| Discriptions | Cat.No. |
|---------------------------|------------|
| UHPLC Guard Column Holder | 5020-08630 |

Replacement Cartridge Ei

| Hardware | Particle Size (µm) | Length (mm) | Phase | ID (mm) | |
|----------|--------------------|-------------|---------|------------|------------|
| | | | | 2.1 | 3 |
| PEEK | 5 | 10 | 100-C18 | 5020-42277 | 5020-42276 |
| | | | 100-SEC | 5020-42330 | - |
| | | | 200-C18 | 5020-42203 | 5020-42202 |
| | | | 300-C18 | 5020-42156 | 5020-42152 |
| | | | 300-C4 | 5020-42158 | 5020-42154 |
| | | | 300-C8 | 5020-42157 | 5020-42153 |
| | | | 300-SEC | 5020-42159 | 5020-42155 |
| | | | HILIC | 5020-42232 | 5020-42231 |

Cartridge Ei Holder

| Discriptions | Cat.No. |
|----------------------------------|------------|
| Ei Holder for Ei Guard Cartridge | 5020-08650 |

Replacement Guard Cartridge E


| Hardware | Particle Size(µm) | Length (mm) | Phase | ID (mm) | | |
|----------|-------------------|-------------|---------|------------|------------|------------|
| | | | | 2.1 | 3 | 4 |
| SS | 3 | 10 | 200-C18 | 5020-42218 | - | - |
| | | | 200-C8 | 5020-42311 | - | - |
| | 5 | 10 | 100-SEC | - | 5020-42331 | - |
| | | | 200-C18 | - | - | 5020-42201 |
| | | | 300-C18 | - | - | 5020-42148 |
| | | | 300-C4 | - | - | 5020-42150 |
| | | | 300-C8 | - | - | 5020-42149 |
| | | | 300-SEC | - | - | 5020-42151 |
| | | | HILIC | - | - | 5020-42230 |

Cartridge E Holder


| Discriptions | Cat.No. |
|-----------------------------------|------------|
| E Holder for 10mm Guard Cartridge | 5020-08500 |

Packed Guard Column


| Hardware | Particle Size(µm) | Length (mm) | Phase | ID (mm) | | |
|----------|-------------------|-------------|---------|------------|------------|------------|
| | | | | 4.6 | 6 | 7.6 |
| SS | 5 | 50 | 300-SEC | - | 5020-42143 | - |
| | | | 100-SEC | 5020-42314 | 5020-42324 | 5020-42329 |



Guard Columns for UHPLC



Cartridge Guard Column Ei (Non-metal Type)



Cartridge Guard Column E



Packed Guard Column

Pre-Column Coupler

Pre-column couplers are joints that are used to connect various guard columns to analytical columns. PCTFE can be used with acids, alkalis, and general organic solvents. In addition, stainless steel can be used under high pressures.

The UHPLC-compatible pre-column coupler is hand-tightened to 50 MPa pressure resistance by means of a special ferrule that is composed of two different materials (i.e., PEEK and metal).

The ferrals are removable and can be used repeatedly, reducing the effect of dead volume when connecting to LC columns with different joint formats.

Specification

Max. operating pressure : 14.7 MPa (Pre-column Coupler UP, Pre-column Coupler W)
 50 MPa (Column Coupler for UHPLC)
 80 MPa (Pre-column Coupler SS)

Tube O.D. : 1/16 inch

| Item | ID (mm) | Length (mm) | Material | Conection | Cat.No. |
|--------------------------|---------|-------------|----------|-----------------------|------------|
| Pre-column Coupler UP | 0.18 | 29.2 | PCTFE | Parker type (UP type) | 6010-49200 |
| | 0.25 | | | | 6010-49201 |
| | 0.50 | | | | 6010-49202 |
| Pre-column Coupler W | 0.25 | 32 | PCTFE | Waters type | 6010-49251 |
| | 0.50 | | | | 6010-49211 |
| Pre-column Coupler SS | 0.10 | 40 | SS | - | 6010-49210 |
| | 0.25 | | | | 6010-49250 |
| Column Coupler for UHPLC | 0.25 | 75 | PEEK, SS | - | 6010-49255 |
| | | 150 | | | 6010-49256 |
| | 0.1 | 75 | | | 6010-49257 |
| | | 150 | | | 6010-49258 |

* Product Pre-column Coupler SS do not fixed with Ferrules, therefore columns with 10-32UNF specification all can be used.



Pre-column Coupler (PCTFE)



Pre-column Coupler SS



Column Coupler for UHPLC

Pre-Column Filter

Specification

Fit in tube O.D. : 1/16 inch

Screw specification : 10-32UNF


Filter pore size : 2 μm

Max. operating pressure : SS : 41.4 MPa (414 bar), PEEK : 34.5 MPa (345 bar)


| Item | Jacket Material | P/N | Qty. (pc) | Cat.No. |
|---|-----------------|-------|-----------|------------|
| Pre-column filter 2 μm | SS | A-315 | 1 | 6010-55100 |
| Replacement pre-column filter 2 μm | - | A-101 | 1 | 6010-55110 |
| PEEK pre-column filter 2 μm | PEEK | A-355 | 1 | 6010-55300 |
| Replacement PEEK pre-column filter 2 μm | - | A-700 | 1 | 6010-55310 |



Pre-column filter



Analysis of Peptides and Proteins




Conditions

Column : ProteoSil 200-C8 (5 μ m, 150 x 4.6 mm I.D.)
 Eluent : A) 0.1% TFA in CH₃CN
 B) 0.1% TFA in H₂O
 A/B = 20/80 – 20 min – 55/45, v/v
 Flow Rate : 1.5 mL/min
 Col. Temp. : 40 °C
 Detection : UV 220 nm
 Injection Vol. : 5 μ L

Analyte

1. Ribonuclease A (0.2 mg/mL)
 2. Insulin (0.2 mg/mL)
 3. Lysozyme (0.2 mg/mL)
 4. BSA (0.2 mg/mL)

Analysis of Protein (ProteoSil 300-C8)




Conditions

Column : ProteoSil 300-C8 (5 μ m, 150 x 4.6 mm I.D.)
 Eluent : A) CH₃CN/0.05 % TFA in H₂O = 80/20, v/v
 B) CH₃CN/0.05 % TFA in H₂O = 10/90, v/v
 A/B = 0/100 – 30 min – 100/0, v/v
 Flow Rate : 1.0 mL/min
 Col. Temp. : 30 °C
 Detection : UV 280 nm
 Injection Vol. : 20 μ L

Analyte

1. DL-Phenylalanine (1.01 mg/mL)
 2. Cytochrome C (0.11 mg/mL)
 3. Lysozyme (0.07 mg/mL)
 4. BSA (0.21 mg/mL)
 5. α -Chymotrypsinogen A (0.08 mg/mL)
 6. Ovalbumin (0.30 mg/mL)

Analysis of Protein (ProteoSil 300-SEC)




Conditions

Column : ProteoSil 300-SEC (5 μ m, 250 x 4.6 mm I.D.)
 Eluent : 0.1 M Na₂HPO₄ (pH 6.9, NaH₂PO₄)
 Flow Rate : 0.5 mL/min
 Col. Temp. : 30 °C
 Detection : UV 220 nm
 Injection Vol. : 20 μ L

Analyte

1. Thyroglobulin (0.25 mg/mL)
 2. BSA (0.25 mg/mL)
 3. Insulin Chain A (0.25 mg/mL)
 4. Oxytocin (0.21 mg/mL)

Analysis of Peptides and Proteins (ProteoSil 300-C18)




Conditions

Column : ProteoSil 300-C18 (5 μ m, 150 x 4.6 mm I.D.)
 Eluent : A) 0.05% TFA in (CH₃CN/H₂O = 90/10, v/v)
 B) 0.05 % H₂O
 A/B = 20/80 – 40 min – 40/60, v/v
 Flow Rate : 1.0 mL/min
 Col. Temp. : 30 °C
 Detection : UV 280 nm
 Injection Vol. : 20 μ L

Analyte

1. Oxytocin (0.05 mg/mL)
 2. Methionine Enkephalin (0.11 mg/mL)
 3. Leucine Enkephalin (0.11 mg/mL)
 4. Angiotensin II (0.05 mg/mL)
 5. Angiotensin I (0.16 mg/mL)
 6. Insulin (0.05 mg/mL)
 7. Insulin Chain B (0.10 mg/mL)

Analysis of Peptides and proteins (ProteoSil 300-C18)




Conditions

Column : ProteoSil 300-C18 (5 μ m, 150 x 4.6 mm I.D.)
 Eluent : A) CH₃CN/0.05 % TFA = 80/20, v/v
 B) CH₃CN/0.05 % TFA = 10/90, v/v
 A/B = 0/100 - 30 min - 100/0 - 10 min - 100/0, v/v
 Flow Rate : 1.0 mL/min
 Col. Temp. : 30 °C
 Detection : UV 280 nm

Analyte

1. DL-Phenylalanine (FW 165)
 2. Cytochrome C (FW 13,000)
 3. Lysozyme (FW 14,000)
 4. BSA (FW 66,000)
 5. α -Chymotrypsinogen A (FW 25,600)
 6. Ovalbumin (FW 45,000)

Analysis of Peptides and proteins (ProteoSil 300-C8)




Conditions

Column : ProteoSil 300-C8 (5 μ m, 150 x 4.6 mm I.D.)
 Eluent : A) 0.1 % TFA in (CH₃CN/0.1 % TFA = 90/10, v/v)
 B) 0.1 % TFA
 A/B = 20/80 - 25 min - 60/40 - 5 min - 60/40, v/v
 Flow Rate : 1.0 mL/min
 Col. Temp. : 30 °C
 Detection : UV 280 nm

Analyte

1. Ribonuclease A (FW 13,700)
 2. Insulin (FW 6,000)
 3. Cytochrome C (FW 13,000)
 4. Lysozyme (FW 14,000)
 5. BSA (FW 66,000)
 6. STI (FW 21,000)
 7. Ovalbumin (FW 45,000)

Analysis of Peptides and proteins (ProteoSil 300-C4)




Conditions

Column : ProteoSil 300-C4 (5µm, 150 x 4.6 mm I.D.)
 Eluent : A) 0.2 % HCOOH in (CH₃CN/H₂O = 90/10, v/v)
 B) 0.2 % HCOOH
 A/B = 20/80 – 20 min – 80/20, v/v
 Flow Rate : 1.0 mL/min
 Col. Temp. : 30 °C
 Detection : UV 280 nm

Analyte

1. Neutrotenin (FW 1673)
 2. Leucin Enkephalin (FW 556)
 3. Cytochrome C (FW 12,000)
 4. Insulin (FW 6,000)
 5. BSA (FW 66,000)
 6. Myoglobin (FW 17,000)
 7. Creatine amidinohydrolase (FW 43,000)
 8. Ovalbumin (FW 45,000)
 9. Creatinine amidohydrolase (FW 170,000)

Analysis of BSA Digests




Conditions

Column : ProteoSil 200-C18 (1.9 µm, 150 x 2.1 mm I.D.)
 Eluent : A) 0.1% TFA in CH₃CN
 B) 0.1% TFA in H₂O
 A/B = 10/90 – 30 min – 50/50 – 0.1 min – 90/10
 – 5 min – 90/10 – 0.1 min – 10/90 – 15 min
 Flow Rate : 0.2 mL/min
 Col. Temp. : 40 °C
 Detection : UV 210 nm
 Injection Vol. : 10 µL

Analyte

Tryptic Digest of BSA (0.5 mg/mL)

Analysis of Oligonucleotides




Conditions

Column : ProteoSil 200-C18 (3 µm, 100 x 3.0 mm I.D.)
 Eluent : A) 5 mM TEAA in H₂O (pH 6.5)/CH₃CN = 80/20, v/v
 B) 5 mM TEAA in H₂O (pH 6.5)
 A/B = 5/95 – (60 min) – 50/50, v/v
 Flow Rate : 0.8 mL/min
 Col. Temp. : 40 °C
 Detection : UV 260 nm
 Injection Vol. : 1 µL
 Sample : Standard

Analyte

1. CATGACGTTCTGATGCT (18 mer, M.W. 5465.61)
 2. CCATGACGTTCTGATGCT (19 mer, M.W. 5754.79)
 3. TCCATGACGTTCTGATGCT (20 mer, M.W. 6058.99)

Analysis of IgG




Conditions

Column Cat. No.: 5020-10818
 Column : MonoSelect RP-mAb (20 x 2.1 mm I.D.)
 Eluent : A) 0.1%TFA in CH₃CN
 B) 0.1%TFA in H₂O
 A/B = 5/95 – 5 min – 90/10 – 0.1 min – 5/95 – 3
 Flow Rate : 0.3 mL/min
 Col. Temp. : 80 °C
 Detection : UV 210 nm
 Injection Vol. : 5 µL
 Sample : Standard

Analyte

1. IgG 0.1 mg/mL

Analysis of Abraxane




Conditions

Column : MonoSelect nPEC (50 x 3.0 mm I.D.)
 Column Cat. No.: 5020-10816
 Eluent : A) CH₃OH
 B) 50 mM Na₂SO₄ in 20 mM Acetate buffer (pH 4.6)
 A/B = 30/70, v/v
 Flow Rate : 1.0 mL/min
 Col. Temp. : Room temperature (approx. 25 °C)
 Detection : UV 260 nm
 Injection Vol. : 10 µL

Analyte

1. Abraxane (Nanoparticle albumin-bound paclitaxel)
 2. Paclitaxel
 * Unknown peak

Analysis of Oligonucleotides



Conditions

Column : ProteoSil HILIC (1.9 µm, 150 x 2.1 mm I.D., Metal-free hardware)
 Eluent : A) CH₃CN
 B) 200 mM HCOONH₄ in H₂O

| Time (min) | A (vol%) | B (vol%) |
|------------|----------|----------|
| 0 | 70 | 30 |
| 45 | 45 | 55 |

Flow Rate : 0.2 mL/min
 Col. Temp. : 40 °C
 Detection : UV 260 nm
 Injection Vol. : 1 µL
 Sample : Standard solution of Single-stranded DNA

Analyte

1. TATGACGTTCTGATGCT (18 mer) 20 nmol/mL
 2. CTATGACGTTCTGATGCT (19 mer) 20 nmol/mL
 3. GCTATGACGTTCTGATGCT (20 mer) 20 nmol/mL



Contact us or your local GL Sciences representative.

<https://www.glsciences.com/contactus/>

Authorized distributor:

<https://www.glsciences.com/company/distributor.html>

For analytical use only. Information subject to change without prior notice.



<https://www.glsciences.com>
E-Mail world@glsc.co.jp

© GL Sciences Inc.
Published in Japan, December 14, 2023
LCBR0303EN

