

Sample Preparation

500 μL sample solution :
First, desalt the peptide solution using MonoSpin C18. Then, dissolve the desalted sample solution with 0.1 % formic acid.

Centrifugation
Speed : 10,000 x g



Centrifuge

30 sec



Centrifuge

30 sec



1. Conditioning


Add 300 μL 0.1 % formic acid

2. Adsorption

Add 500 μL peptide sample solution

3. Elution

Sample solution in 0.1 % formic acid



Always replace and attach a new recovery tube whenever adding a new elution buffer

Details of each elution buffer

① 25 mM HCOONH ₄	200 μL
② 50 mM HCOONH ₄	200 μL
③ 100 mM HCOONH ₄	200 μL
④ 500 mM HCOONH ₄	200 μL
⑤ 1 M HCOONH ₄	200 μL

* Each elution buffer contains 10 % acetonitrile

Conditions

Column : Inertsil ODS-3 (3 μm , 2.1 x 150 mm)
Eluent : A) H₂O (0.1 % HCOOH)
 B) Acetonitrile (0.1 % HCOOH)
 A/B = 90/10 - 20 min - 50/50
Detection : UV 210 nm
Flow Rate : 0.2 mL/min
Col. Temp. : 40 °C
Injection Vol. : 2 μL