

**anthos**

**A3**



**Anthos  
Classe  
A3Plus**

CONTINENTAL  
INTERNATIONAL

**improve**

## CLASSE A3 PLUS CONTINENTAL INTERNATIONAL

**im**

**prove.**



### **improve potential.**

Classe A3 Plus : Anthos world improves your potential by adding the extra safety, reliability, flexibility and functionality that are key to meeting every clinical need.

### **improve quality.**

Outstanding construction quality, cutting-edge technology, excellent operating ergonomics, patient comfort. Classe A3 Plus, a clever investment to secure advantages well beyond the usual standards.

### **improve performance.**

Classe A3 Plus offers a wide range of specific operational features that raise the bar of the dental unit's performance profile. The unit's features can be customised for maximised results.

improve potential.  
improve quality. improve  
performance.

## CLASSE A3 PLUS

The Digit control panel layout, complete with X-ray viewer, is easy to read and includes several features. Access to the settings of each individual instrument, patient chair and integrated devices is fast, user-friendly and can be viewed on the alphanumeric display unit. Its design maximises operational ergonomics, making the operator's work easier.



CE  
0051

ISO 9001:2008 EN ISO 13485:2012

### Simple, intuitive operation

Touch screen display.

The control panel surface, protected by capacitive glass, can be easily disinfected. Quick adjustment of Speed/Power with numerical indication & other signals (bottle, light etc.) Clean button is also available




### Maximised Adaptability

Continental :

- ✦ The redesigned dentist's module and arms system allows the operator to work in any operating position. Classe A3 Plus responds ergonomically to every clinical need
- ✦ Up to 6 instrument arms
- ✦ Each instrument has retraction valve system to avoid cross contamination



New instrument lever design - with shorter vertical footprint to reduce interference with the operating light and ensure increased instrument outreach.



ISO-JOINT Ergonomics / Trendelenburg Position : New ISO-JOINT geometry, compensated backrest-seat movement maximises comfort and minimises sliding of the patient's head. Easy to maintain Trendelenburg position in case of any emergency. Easy maintenance of Cradle position enhance comfort of patient.

Vertical travel and lifting capacity. A long vertical travel (37-80 cm) makes it easier to provide treatment in comfort and allows easy patient chair access for the elderly. The patient chair is certified to lift patients weighing up to 190 kg.



Operational ergonomics and patient comfort are the natural outcome of the unit's design.





**Improved Productivity**

Micromotors (Optic) :

The supplied model i-MMr (3.3 Ncm) is also available. Speeds can set between 100 and 40,000 rpm.

Able to cover every modern dentistry need, light weight, balanced fiber optic lighting.

**Scalers (Optic) :**

With LED handpieces compatible with the best tips on the market. Effective & Safe, excellent for supra & subgingival application.



**Condensate drain at chair front :**

For absorption of extra moisture in the compressed air & enhance life of the chair.

**Syringes:**

Ergonomically-shaped 3 way syringe are available in delivery unit as well as assistant arm.



**TFT Clamp with Multimedia Cable :**

Provides easy installation of monitor for flexible dental operation

**Headrest :**

In addition to an adjustable 2-axis version with mechanical lock.



## Protected Environment

Safe clinical environment for the patients, the operating team and the dentist, thanks to a wide range of integrated, optional hygiene systems. Automatic devices that prevent contamination and design features that minimise the risk of cross-contamination.



### Filter :

Large double-chamber filter ensuring greater suction system efficiency

### Dentist's Module Handle :

The Continental dentist's module handle can be removed and autoclaved



### Adjustable Knob in the Delivery Unit Arm

For easy adjusting while performing dental procedure & prevent any displacement of the unit.

### Chair Safety

Safety feature at chair base, backrest & assistant arm to avoid any unwanted mishap.



### Cuspidor Bowl Unit

Ceramic spittoon, it can be completely disassembled for effective sanitization. Removable & disinfectable cup filler & cup holder units.

## Complete Solution



### Removable Support :

On the Continental module, the instrument support can be removed and disinfected to avoid cross contamination.

### Upholstery

Easy-to-sanitise seamless, durable upholstery. It is softer & more comfortable. Greater resistance to tears & dirt.



### Foot Control

Foot controls with multiple functions. These allow activation of Zero position / Chip air & chip water / Joystick for dental chair movement / Instrument speed & power adjustment i.e. pressure operated foot control / rinsing position.

### Stop Vacuum

Device integrated in the patient chair base: when pressed it stops suction without the user having to replace the cannulae in the holders.



### Venus Plus L-LED

With hermetically sealed front screen, enables on-off control also via an infrared sensor and features adjustable light intensity from 3,000 to 50,000 Lux with 5000K colour temperature. Power adsorption- low

### Assistant Arm

3 Slot Assistant Arm : low & high vacuum suction. Control Panel on Assistant side for easy operation.





WHERE **TRUST** COMES **FIRST**



**Warranty Ke Saath Bhi  
Warranty Ke Baad Bhi**

**For any Sales & Support**

 98211-84441, 84483-98852

 [sales@unicorndenmart.com](mailto:sales@unicorndenmart.com)

**Customer Care**

Register Your Complaint with A Missed Call: **7065-558-558**

 011-45551250, 51, 55  9717012081, 95994 48518

 [care@unicorndenmart.com](mailto:care@unicorndenmart.com) | [care1@unicorndenmart.com](mailto:care1@unicorndenmart.com)

Connect With



Unicorn DenMart



Corporate Office: 3 Local Shopping Centre, MOR Land, Near J Block DDA Market,  
New Rajender Nagar, New Delhi – 110060



\* Accessories shown may not be a part of standard equipment. Specifications are subject to change without notice.

Ver 5/2026