



HUNAN BEYOND MEDICAL TECHNOLOGY CO., LTD


Addr: Beyond Zone,Lijiacun Rd,Xueshi Street,Yuelu District, Changsha, 410208,Changsha,China


Contact person: Sandy Chu


Mobile:+86-17352837131



E-mail: sandy@csbeyond.com

Quotation

Model no.	Specification	Picture	FOB Price(USD)		
			5-19 pcs	20-49 pcs	>=50 pcs
Endo motor	<ol style="list-style-type: none">1. Dimensions: 140x220x92mm2. Weight: 460 grams3. Rated voltage and frequency: Input: ~220V, 50Hz; Output: DC 5V 1.0A Output speed range: 150-625 r/min4. Output torque range: 0.8-3.2 N.cm5. Mobile phone speed ratio: 16:16. Internal battery: Rechargeable battery (lithium battery) Specification: 3.7V, 1050mAh		110	95	85

<p>LED Curing Light</p>	<p>1. Features Can cure various resin materials, especially suitable for orthodontic brackets, porcelain veneers and core post adhesives.</p> <p>2. Wireless charging, large-capacity lithium battery, after fully charged, irradiate with the maximum radiation intensity, can be used continuously for more than 500 times (each time is less than 10 seconds)</p> <p>3. The battery can be replaced, which can effectively extend the service life of the machine</p> <p>4. Time setting 5s, 10s, 15s, 20s, 25s, 30s, 35s, 40s</p> <p>5. Mode selection Three working modes of continuous irradiation, pulse and soft start are optional</p> <p>6. The highest light intensity can reach 1600mw/cm²</p> <p>7. Spectral range: 420~480nm</p>		<p>25</p>	<p>20</p>	<p>15</p>
-------------------------	--	--	-----------	-----------	-----------

<p>GP Obturation system</p>	<ol style="list-style-type: none"> 1. Two-in-one filling system, convenient, fast and accurate filling and hot pressing; equipped with four heating needles of different specifications and two glue injection needles of different lengths. 2. Exquisite design, simple operation and ergonomic design; 3. Simple technology and easy to learn; 4. The heating needle of the hot pressing pen is heated to the set temperature within 1 second, and it is operated when it is pressed and heated; 5. The heating pen is heated by 360° ring buttons, with a lightweight body, comfortable holding and convenient operation. 6. The temperature of the filling gun is displayed in real time, and the operation is simple and clear; 7. Wireless charging method, non-contact contact, ready to use with charging, convenient and fast; full charge, the filling gun can work continuously for more than 60 minutes at 220 degrees Celsius, and the hot press pen can work continuously for more than 800 times (less than 10 times each second) 8. The filling gun has an automatic energy-saving mode to effectively extend the battery working time. If it is not working for 1 minute, it will automatically enter the energy-saving mode. When the doctor grabs the gun body again to work, it will automatically quickly heat to the set temperature. 		<p>300</p>	<p>260(>10pcs)</p>	<p>240(>20pcs)</p>
-----------------------------	---	---	------------	-----------------------	-----------------------

<p style="text-align: center;">Endo files</p>	<ol style="list-style-type: none"> 1. Brand new imported nickel-titanium material, the file body has the function of restoring the material memory, providing visual warning for the safety of the file 2. The unique heat treatment technology can treat the file body to any hardness; the blue file body is softer and tougher, and the yellow file body is harder and stronger in cutting force. 3. The wise variable cone design preserves the dentin to the greatest extent during the effective root canal enlargement process, achieving minimally invasive purposes. 4. Reasonable thread design, effectively avoid the risk of file breaking; superb manufacturing technology, perfect process flow, super cutting ability. 		5	4(>500pcs)	3 (>1000pcs)
<p style="text-align: center;">Orthodontic pliers</p>	<ol style="list-style-type: none"> 1. New anti-rust treatment、anti-discoloration、anti-rust 2. Medical grade stainless steel, no plating and polished surface design, comfortable and durable 3. The beak of the pliers is inlaid with imported alloy materials, precision numerical control technology, anti-dropping, wear-resistant. 4. 15 fine manufacturing processes and strict quality inspection to ensure the high quality of each orthodontic forceps. 		See detailed price list		

1. Delivery time

Depends on quantity (Regular: 3 working days)

2. Payment Form

100% TT in advance

3. Warranty

1) Two years warranty since delivery.

2) Repair in factory: Beyond Medical pays the repair fee & delivery, distributor pay the fee of sending back.

4. Price

The FOB price is for the Distributors, based on Real Time Exchange Rates, including the domestic ordinary Freight

5. Packaging

Each in carton box, the box is suitable for sea/air transportation.