



Shockwave



Exclusive to OrthoCanada



BTL-6000 Radial Shockwave

BTL-6000 ELITE Shockwave

Offering the most powerful radial shockwave device on the market with an intensity of up to 6 bars and ultralight Intelligent MAGNUM applicator, the BTL-6000 Elite shockwave unit is a clinical solution for sub-acute and chronic musculotendinous conditions. The intensity of 6 bars will deliver successful treatment of deeper tissues while the ultralight Intelligent MAGNUM applicator accommodates parameter adjustments seamlessly.

MDALL 106498

- 8.4" colour touchscreen
- 22 preset protocols
- Includes Magnum applicator
- Intensity up to 6 bars
- Frequency up to 22 Hz



BTL-RSWT-ELITE BTL-6000 ELITE Shockwave

BTL-6000
ELITE Shockwave

Mechanism of Action

Radial shockwave is an acoustic wave that carries high energy to myoskeletal tissue. It causes mechanical stimulation that promotes tissue regeneration and repair in injured tendons and other soft tissue. It is the first-choice treatment for chronic conditions such as superficial tendinopathies, calcifications, and persistent pain of musculotendinous origin.



Shockwave Training Canada

For a complete schedule of courses and workshops please visit www.shockwavetraining.ca



Intelligent Magnum Applicator

Spin to fine-tune therapy parameters

6 Bars & 22Hz

Shock-absorbing cushioning eliminates reverse shocks

Smart indicator of therapy parameters

Ultralight & ergonomically-shaped

Redesigned projectile construction for extended lifespan



BTL-6000 PRO Shockwave

- 5.7" colour touchscreen
- 22 preset protocols
- Includes Magnum applicator
- Intensity up to 4 bars
- Frequency up to 15 Hz



BTL-6000 EASY Shockwave

- User friendly LED interface
- 7 preset protocols
- Intensity up to 4 bars
- Frequency up to 15 Hz

Common Indications



Patellar Tendinitis



RC Tendinitis, Bursitis & Calcification



Lateral Epicondylitis



Trigger Point Release



Hamstring Strain



BTL-6000 Focus Shockwave

The BTL-6000 Focus Shockwave device is a clinically proven and effective non-invasive alternative to surgery. Ideal for fast pain relief and the restoration of mobility, the therapy enhances tissue regeneration, eliminates calcification, relaxes muscles and trigger points, and stimulates the healing of non-union fractures. Its technology offers non-decreasing shock intensity over its lifespan, supported by constant shock intensity monitoring. Optimal Focal Zone covers any treated area easily and the high Energy Flux Density enables treatment of diverse tissues. It has the highest frequency range on the market for a wide spectrum of indications.

MDALL 102793

BTL-FSWT

BTL-6000 Focus Shockwave



BTL-6000 Focus Shockwave with Cart

Medical Effects

Acoustic waves with high energy peak used in Shockwave Therapy interact with tissue causing overall beneficial medical effects – accelerated tissue repair and cell growth, analgesia, and restored mobility. Each of the following effects happen simultaneously and are used to treat chronic, sub-acute and acute (advanced users only) conditions.

- New blood vessel formation
- Reversal of chronic inflammation
- Stimulation of collagen production
- Dissolution of calcified fibroblasts
- Dispersion of pain mediator substance P
- Release of trigger points



Common Indications



Tendinitis & Bursitis



Non-Union Fracture



Calcification



Trigger Point Release



Avascular Necrosis

Most Advanced Features in the Industry

Patented electroacoustic technology offering the highest lifespan of 3 million shocks

Intelligent applicator with continuous monitoring of shock intensity

Adjust parameters directly from the applicator touchscreen display

Easy maintenance & liquid-free circuit

Switch & go applicator replacement service



Select the depth of penetration of up to 65 mm with 3 different coupling pads



Your OrthoCanada Sales Consultants



Zach El-Shafei
zach.elshafei@orthocanada.com
1-800-561-0310 ext 8000



Karen Wright
karen.wright@orthocanada.com
1-800-561-0310 ext 1719



Razi Khaddage
razi.khaddage@orthocanada.com
1-800-561-0310 ext 1704



Jonathan Hamel
jonathan.hamel@orthocanada.com
1-800-561-0310 ext 2020

