


InertSearch™ for GC

InertCap® Applications

Volatile organic compounds

Data No. GA005-0000



System : GC/MS P&T/GC
注) 高速分析用 HTP セットが必要になります。

Column : AQUATIC-H

Injection : P&T

Detection : MS

Solutes : 1. 1,1-Dichloroethene
2. Dichloromethane
3. t-1,2-Dichloroethene
4. cis-1,2-Dichloroethene
5. Chloroform
6. 1,1,1-Trichloroethane
7. Carbone Tetrachloride
8. 1,2-Dichloroethane
9. Benzene
10. Trichloroethene
11. 1,2-Dichloropropane
12. Bromodichloromethane
13. cis-1,3-Dichloro-1-propene
14. Toluene
15. t-1,3-Dichloro-1-propene
16. 1,1,2-Trichloroethane
17. Tetrachloroethene
18. Dibromochloromethane
19. m-Xylene
20. p-Xylene
21. o-Xylene
22. Bromoform
23. p-Dichlorobenzene
I.S. p-Bromofluorobenzene