

THE PERFECT **ENDO**



Wave Q

Sonic Irrigation Activator

Wave Q is a specially designed Endo Activator for dentists. With the working tip, the flushing solution is activated through the working tip's vibration to clean the root canal.

13,000

vibrations per minute

62g

holding weight

90 days

battery duration

20mm

holding size



Wave Q

	SIZE	176 X 22 X 22mm
	WEIGHT	60±2G
	BATTERY	3.7V / 800mAh lithium battery
	OUTPUT	5V DC, 2A
	INPUT	AC 100-240~50/60HZ
	MATERIALS	Contra angle: plastic and silicone Working tip: POM plastic

Features

Easy to use

3 different sonic cleaning frequencies.
Easy to install and disassemble in one step.

Avoid breakage


Polymer material(POM) material working tips
avoid breakage during cleaning.

Ergonomic design

Only 62g and the grip part diameter is 20mm
It is the most comfortable size for holding.

Safe and hygienic

Steam sterilizable contra angle avoids cross-infection
Sterilization condition: 135°C, 4 minutes, 0.22Mpa



STREAM

Traditional anesthetic procedure is manual controlled and the extrusion speed is unstable, which made the poor absorption effect and the painful feeling of patient.

Stream X is an assistor booster which is specially designed in length and weight, the length is small and the center of gravity of the product is in the front. It relieves patients' pain and doctors' fatigue.

MODE

TIME

Fast Mode

1min

Slow Mode

2min

PDL Mode

4min

Custom Mode

1-4min

Length

166
mm

Weight

72
g

Reflow Time

30
s

Battery Duration

30
days

ASSISTOR
BOOSTER



STREAM X

FEATURES

★ Custom Injection Setting

Stream X is equipped with 4 modes of injection time which varies from 1 to 4 minutes.

★ Short Reflow Setting

The reflow time is just 30 seconds. It's really convenient for the doctor to start the next operation.

★ Comfortable to Hold


Ergonomic & comfortable designed, weight 72g, holding size 20mm, the length is small and the center of gravity is in the front.

★ Durable and Efficient

The 350mAh lithium battery can be continuous used for 30 days(6 root canals per day).

SPECIFICATION

SIZE	166x21x33mm
WEIGHT	73±2g
BATTERY	3.7V/350mAh lithium battery
BATTERY DURATION	30 days
OUTPUT POWER	5V DC, 2A
INPUT POWER	100-240~50/60Hz
CHARGING TIME	2h(Recommended to charge for more than 3 hours for the first time.)



GP

Gutta Percha
Cutting Pen

The GP Flow is specially designed for dentists to soften or cut the gutta-percha during RCT by heating the working tip.



FLOW

GP FLOW

Gutta Percha Cutting Pen

Features

Efficient

Rapid heating in 0.2s at 4 levels temperature: 150°C, 180°C, 200°C, 230°C

Easy to use

The three-button operation logic is simple, making it easy to use.

Durable

Equipped with 2600mAh large battery and can be continuous used for 30 days

Comfortable

Ergonomic design: comfortable holding weight 95g and size 20mm, and the working tip is 360° rotatable.

Information

SIZE	174 x 25 x 25mm
WORKING TIP SIZE	40/.04; 45/.04; 50/.04; 55/.08
HEATING TEMPERATURE	150°C, 180°C, 200°C, 230°C
WEIGHT	95±2g
BATTERY	3.7V / 2600mAh lithium battery
OUTPUT POWER SUPPLY	5V DC, 2A
INPUT POWER SUPPLY	AC 100-240~50/60Hz
STERILIZATION CONDITIONS	135°C, 4minutes, 0.22Mpa





WHERE **TRUST** COMES FIRST



**Warranty Ke Saath Bhi
Warranty Ke Baad Bhi**

For any Sales & Support

 9599281852, 8448398852

 sales@unicornndenmart.com



SCAN TO VISIT OUR WEBSITE



Corporate Office: 3 Local Shopping Centre, MOR Land, Near J Block DDA Market,
New Rajender Nagar, New Delhi – 110060

    /unicornndenmart.com

* Accessories shown may not be a part of standard equipment. Specifications are subject to change without notice.

Ver 03/2025