


InertSearch™ for GC

InertCap® Applications

VOC & Amines

Data No. GA065-0000



System : GC/FID
Column : InertCap for Amines
0.32mmI.D. × 30m
Column Temp. : 50°C
Carrier Gas : He 70kPa
Injection : Split 1 : 50
Detection : FID

Solutes : 1. Methanol
2. Dimethylamine (DMA)
3. Trimethylamine (TMA)
4. Ethanol
5. iso-Propanol
6. n-Pentane
7. n-Propanol
8. n-Hexane