KTape Always a Step Ahead

The Original



Exclusive to OrthoCanada



Always a Step Ahead

Recommended by the K-Taping Academy Canada, K-Tape is an innovative product that has gained global recognition for its exceptional quality and performance. The Kumbrink-developed K-Tape combines the best natural, raw materials with the most advanced manufacturing process available. Its impact on generating positive clinical outcomes has been profound, introducing a wide array of new possibilities for therapeutic interventions across many healthcare providers.

Cotton & Dyes

K-Tape uses high-quality cotton from a single supplier ensuring uniform thread density and knot-free manufacturing. To maintain consistent elasticity throughout its use, the tape incorporates premium elastane thread while dyes used are supplied by Dystar, known for their safety and adherence to the OEKO-TEX Standard 100 label.

Backing Paper

K-Tape is pre-stretched to the optimal level of 10% during its application to the backing paper. This ensures that it has the ideal amount of pre-stretch required for muscle technique taping.

Physiobond Adhesive

K-Tape features professional-grade Physiobond adhesive, exclusively used for K-Tape. The adhesive undergoes a meticulous manufacturing process, creating long polymer chains for effective skin adhesion. A complex temperature-time procedure removes any skin-irritating residual monomers, resulting in a top-notch tape that sets a new industry standard.

K-Taping

K-Taping is a non-pharmaceutical therapy technique that is effective for an extremely wide range of clinical conditions. K-Tape's high-quality fabric is water resistant and breathable, and designed to provide longer wear and comfort. K-Tape is pre-stretched and is ideal for muscle technique taping. Professional training is essential for effective application of K-Tape therapy. For more information visit www.k-taping.ca.







K-Tape – Box of 4 Rolls

OC-4343	А	Blue	2" x 16.4' each
OC-4342	в	Red	2" x 16.4' each
OC-4301	С	Beige	2" x 16.4' each
OC-4341	D	Black	2" x 16.4' each

K-Tape - Mixed Colours

OC-4344 E	4 Rolls	2" x 16.4' each

K-Tape – Bulk Roll

OC-4373	F	Blue	2" x 72'
OC-4372	G	Red	2" x 72'
OC-4365	н	Beige	2" x 72'
OC-4371	I	Black	2" x 72'

K-Tape Pure Cotton

OC-4364	J	Bulk Roll	2" x 72'
OC-4331	К	Box of 4 Rolls	2" x 16.4' each

K-Tape My Skin - Box of 4 Rolls

OC-4301	L	Beige	2" x 16.4' each
OC-4302	М	Light Brown	2" x 16.4' each
OC-4303	Ν	Medium Brown	2" x 16.4' each
OC-4304	0	Dark Brown	2" x 16.4' each
OC-4305	Р	Very Dark Brown	2" x 16.4' each

K-Tape My Skin – Mixed Colours

OC-4300	Q	Box of 5 Rolls	2" x 16.4' each

K-Tape My Skin – Bulk Roll

OC-4365	R	Beige	2" x 72'
OC-4366	S	Light Brown	2" x 72'
OC-4367	т	Medium Brown	2" x 72'
OC-4368	U	Dark Brown	2" x 72'
OC-4369	V	Very Dark Brown	2" x 72'



K-Tape Team Colours

Display your team spirit by proudly wearing your team colours. Customize your K-Tape selection across 9 different colour options.



K-Tape Starter Kits

Save with K-Tape starter kits. They Include a box of (A) 4 mixed colour K-Tape rolls, a pair of either (D & E) 21cm K-scissors or (F) 16cm bandage scissors, an (c) 85ml tube of Pre K Gel, and a (g) package of Crosstape mixed sizes.

OC-4358	Kit #1	w/21 cm K-Scissors Right Hander
OC-4357	Kit #1	w/21 cm K-Scissors Left Hander
OC-4359	Kit #2	w/16 cm Bandage Scissors Right Hande



Kumbrink K-Taping 2nd Edition

This fully illustrated guide to K-Taping by Birgit Kumbrink, BSc PT is an industry standard and best seller. Includes 450 photos, detailing K-Taping techniques for 50+ conditions across 200 pages. It describes the efficacy of K-Taping in healing, pain relief, and biomechanics. An essential for courses and therapists.

Pre K Gel

Used for skin preparation prior to the application of K-Tape. Its alcohol based formula cleans and disinfects the skin, while the antiperspirant properties of aloe vera and green tea improve adhesion and duration of wear.

OC-4374 C 85 ml Tube

K-Tape Scissors

Blades are teflon-coated to prevent the acrylic adhesive infiltrating the blade's metallic pores and to retain long lasting sharpness. Manufactured in Japan.

OC-4353	D	21 cm K-Scissors Right Hander
OC-4352		21 cm K-Scissors Left Hander
OC-4354	F	16 cm Bandage Scissors Right Hander
OC-4361	Е	13 cm K-Scissors Pink with Band Clip

CROSSTAPE[®]


Therapy for Trigger & Pain Points


The Kumbrink Crosstape patch is used on trigger points, acupuncture points, tense muscles and sore joints. It can effectively be used to stimulate selected acupuncture points, enabling an acupuncture-like treatment without the use of needles.

Many Sizes. Many Applications.


The Crosstape patch is available in sizes S, M, L and XL. Size M is the most commonly used. Size S is ideal for facial applications. Sizes L and XL can be used effectively on larger areas.

OC-4351-S	А	Small	400 Patches
OC-4351-M	в	Medium	180 Patches
OC-4351-L	С	Large	120 Patches
OC-4351-XL	D	X-Large	40 Patches
OC-4351-MIX		Mix Box	20 Sml, 27 Med, 6 Lrg & 2 XL Patches





Actual sizes illustrated



Common Crosstape Applications

Temporal Headache









Epicondylitis (Tennis Elbow)

0



