

# Light Cure Units



# Introducing


The All New

## Features

- Special lens condenser design
- Broad spectrum blue-violet light with 12W high power LED light source
- 360° rotation curing head
- 2 working modes:  
Strong light mode: 1s,3s,5s  
Standard mode: 10s,20s,30s
- 7 Magnetic special lens
- 6 Month warranty



Seven Magnetic Special Lens



## LIGHT CURE UNIT

### O-Light

**1** Sec  
Curing



- ✦ 1 Sec Curing of Composite Resin
- ✦ High power LED Light with imported wick from America
- ✦ High-intensity light output with strong penetration
- ✦ 360° rotatable probe
- ✦ High power Visible Blue LED light
- ✦ Wavelength: 420 - 490 nm
- ✦ Modes: high mode: 1 sec, 3 sec
- ✦ Standard Mode: 5 sec, 10 sec, 15 sec, 20 sec
- ✦ High mode for fast curing of post & core, orthodontic brackets, veneers in intensity 2300 mw/cm<sup>2</sup>
- ✦ Standard mode for curing in restoration cases in intensity 1200 mw/cm<sup>2</sup>

## *LIGHT CURE UNIT* **TULIP DIGITAL LED**




- ✦ Material of shell handle: Aluminium
- ✦ Size: Diameter 23 X 160mm
- ✦ Mode: Three kinds of working mode
- ✦ Working Volt: 4.2~4.5 volts
- ✦ Electric current: 0.7 amperes
- ✦ Battery Type: Lithium-Ion
- ✦ Battery capacity: 1400Mah
- ✦ Wave Length: 420-480nm
- ✦ Luminous Intensit: 1200mW/cm<sup>2</sup> 5W

## LIGHT CURE UNIT


### O-Light Plus

- ✦ Upgraded wide spectrum  
385nm-515nm
- ✦ Poly wave Blue-violet light, it can cure all types of resins
- ✦ New style resins material contains TPO (photo initiator), can also be cured
- ✦ 360° curing - reach every angle
- ✦ Modes: 1> High Mode - Suitable for Orthodontic Brackets, Porcelain Veneers, Post and 2> Standard Mode - Suitable for General Dental Restorative Procedures
- ✦ Light Intensity: High Mode: 2300~2500 mw/cm<sup>2</sup>  
Standard Mode: 1000~1200 mw/cm<sup>2</sup>
- ✦ Time Settings: High Mode: 1 Sec, 3 Sec Standard Mode: 5 Sec, 10 Sec, 15 Sec, 20 Sec


**1** Sec  
**Curing**




Metal head prevents  
Anti-shock & cracks



Upgraded focused  
light output



O-Light plus Spectral Curve  
& Resin Spectral Absorptivity



Light hood prevents  
eyes from Blue-ray  
injury



**WHERE TRUST COMES FIRST**



*Partnering with Dentists for over 2 decades - together creating a world full of healthy smiles.*

**For any Sales & Support**

 9821184441, 9599281852, 8448398852

 [sales@unicorndenmart.com](mailto:sales@unicorndenmart.com)



SCAN TO VISIT OUR WEBSITE



Corporate Office: 3 Local Shopping Centre, MOR Land, Near J Block  
DDA Market, New Rajender Nagar, New Delhi – 110060

 [www.unicorndenmart.com](http://www.unicorndenmart.com)

\* Accessories shown may not be a part of standard equipment. Specifications are subject to change without notice.

Ver 12/2022