


Dicarboxylic acids

Data No. LA110-0000



Column : Inertsil ODS-2
5µm
150 × 4.6 mmI.D.

Column Cat. No. : 5020-01101

Eluent : A:CH₃CN/0.1% H₃PO₄=2/98
B:CH₃CN/0.1% H₃PO₄=45/55
A/B 100/0?(25min)?0/100

Flow rate : 1.0 mL/min.

Col. Temp. : 40°C

Detection : UV 210 nm

Injection Vol. :

Sample :

Solutes : 1. Malonic acid
2. Succinic acid
3. Glutaric acid
4. Adipic acid
5. Pimelic acid
6. Azelaic acid
7. Sebacic acid