



- 03 Colposcopes
  - **1** LEEP Precision System
    - 05 LEEP Electrodes
      - **06** LEEP System Accessories & Instruments
        - Milex® Pessaries for Prolapse & Incontinence Endocell® Disposable Endometrial Cell Sampler





## LEISEGANG OptiK® Model 2

Optimized Digital Imagery Paired with Exceptional Optics

LEISEGANG celebrates more than 60 years as a world leader in colposcope design and production, recognized for German precision, clarity and durability. The Colposcopes of OptiK Model 2 series are Leisegang's popular Digital/Photo colposcopes with precision fine adjustments in horizontal and vertical direction, high performance optics and LED lighting. Specially selected, hand-polished and hand-set glass optics are used in all LEISEGANG OptiK colposcopes. The optical brilliance achieved through our precision manufacturing process is unparalleled.

#### **FEATURES**

- OptiK scope head with convergent optical beams
- Par-focal, 3-step magnification (3.75x, 7.5x, 15x)
- Bright, white, energy-efficient Light-Emitting Diode (LED) light, stable color spectrum (similar to daylight) provides highly improved contrast and better clarity
- LED life > 20,000 hours or 20 years1 for lasting performance
- Independently adjustable oculars and diopters
- Measurement reticules for quick lesion size assessment
- · Gross/fine focus and height adjustment
- EMR compatibility captures jpgs to add to EMR system







# LEEP PRECISION™ Integrated System

Office Diagnostics and Procedures

The preferred choice for functionality, performance & precision for Loop Electrosurgical Excision Procedures.

The LEEP PRECISION™ Integrated System has been designed to meet the unique requirements of the clinic and practice-based facility. Fully integrated and mobile, the LEEP PRECISION Integrated System combines three components – LEEP PRECISION Cart, LEEP PRECISION Generator, and the LEEP PRECISION Smoke Evacuator – for maximum efficiency.

#### **FEATURES**

- Hand or foot activation for operating precision
- Smoke evacuator eliminates odor at 99.999 percent efficiency level
- Fully integrated LEEP station with storage
- Safe and easy to operate
- Heavy-duty casters facilitate mobility
- Elegantly designed for the modern procedure room








## **LEEP Electrodes**

Tungsten wire electrodes provide exceptional quality & performance. Includes various sizes of loops, square, balls & needle tips specifically designed for loop electrosurgical excision procedures.



#### **Large-Radius Loop Electrodes** ≥ 2 cm Wide

NUMBER	DESCRIPTION
R2510	2.5 cm W x 1.0 cm D • 12 cm shaft
R2020	2.0 cm W x 2.0 cm D • 12 cm shaft
R2015	2.0 cm W x 1.5 cm D • 12 cm shaft
R2010	2.0 cm W x 1.0 cm D • 12 cm shaft
R2008	2.0 cm W x 0.8 cm D • 12 cm shaft

#### **Medium-Radius Loop Electrodes** <2 cm Wide

NUMBER	DESCRIPTION
R1515	1.5 cm W x 1.5 cm D • 12 cm shaft
R1507	1.5 cm W x 0.7 cm D • 12 cm shaft
R1505	1.5 cm W x 0.5 cm D • 12 cm shaft
R1010	1.0 cm W x 1.0 cm D • 12 cm shaft
R1007	1.0 cm W x 0.7 cm D • 12 cm shaft

# **Square Electrodes**

NUMBER	DESCRIPTION		
S1010	1.0 cm W x 1.0 cm D 12 cm shaft		
S1007	1.0 cm W x 0.7 cm D 5.5 cm shaft 15° bend		
\$1004	1.0 cm W x 0.4 cm D 5.5 cm shaft 15° bend		
V1004	1.0 cm W x 0.4 cm D 8 cm shaft		

### **Needle Electrodes**

N1200

V1004

NUMBER	DESCRIPTION
N1200	0.8 mm D x 16 mm 12 cm shaft
N0550	0.8 mm D x 16 mm 5.5 cm shaft
N2512	0.25 mm D x 8 mm 12 cm shaft
N2555	0.25 mm D x 5 mm 5.5 cm shaft

N2512

## **Ball Electrodes**

N2555

ER	DESCRIPTION		NUMBER	DESCRIPTION
10	0.8 mm D x 16 mm 12 cm shaft		B0512	5 mm D 12 cm shaft
0	0.8 mm D x 16 mm 5.5 cm shaft		B0555	5 mm D 5.5 cm shaft
2	0.25 mm D x 8 mm 12 cm shaft		B0312	3 mm D 12 cm shaft
5	0.25 mm D x 5 mm			







## **LEEP System Accessories & Instruments**

Includes all consumables needed for safe & effective electrosurgical conductance as well as for efficient smoke evacuation.

#### Handpiece, Single-Patient Use

The sterile, single-patient use handpiece works with any of our disposable electrodes.
 The wand is operated with a footswitch.

#### **Return Pad, Single-Patient Use**

Non-sterile patient return pads.

#### **Tubing Set, Single-Patient Use**

• Each set includes 6 feet of sterile disposable tubing and a DSE (Disposable Smoke Evacuator) tube. This set completes the connection between the reducer and the specula.

### Euro-Med® LEEP Instruments

A range of nonconductive, fully autoclavable specula offered in a variety of styles to facilitate cervical tissue excision & removal.

#### **Graves LEEP Specula**

The Graves LEEP Specula is designed for electrosurgical applications. Also features integrated smoke tube. Available in three sizes.

Small W16 mm L84 mm Medium W22 mm L114 mm Large W26 mm L121 mm

### **Pederson LEEP Specula**

Designed specifically for electrosurgical applications, the narrow blades of the Pederson are more comfortable for smaller patients or those with a narrow vaginal canal.

 Small
 W13 mm L81 mm

 Medium
 W22 mm L114 mm

 Large
 W25 mm L120 mm









## Milex® Pessaries for Prolapse & Incontinence

Provides a nonsurgical, noninvasive solution for assessing or resolving pelvic organ prolapse & stress incontinence. These silicone devices are available in a wide variety of configurations to provide an optimal fit to meet the specific gynecological & urinary issue.



**A. Hodge**For stress urinary incontinence for a woman who has a very small vaginal introitus.



**B. Incontinence Dish**For stress urinary incontinence with mild prolapse.



C. Incontinence Dish w/Support
For stress urinary incontinence with mild prolapse and a mild cystocele.



D. Gehrung with Knob/Folding For cystocele with or without mild rectocele complicated by stress urinary incontinence.



E. Gehrung/Folding Support for cystocele and rectocele. It is also more adaptable in cases

of procidentia where the uterus tends to herniate when other pessaries are used.





# **F. Cube** Support for 3rd degree

prolapse/procidentia. Should be removed nightly.

# Endocell® Disposable Endometrial Cell Sampler

The Endometrial Cell Sampler is a single-use, sterile, disposable plastic device designed to extract sample tissue from the uterine mucosal lining, or to extract a sample of uterine menstrual content, for histological examination.





# The **Future of Medical Care** in the Caribbean.



\*TRINIDAD & TOBAGO – HEAD OFFICE Corner Churchill Roosevelt H/way & Tissue Drive, Trincity, Trinidad & Tobago Tel: 1 (868) 640 2482 T&T Toll Free 800 4225 Email: enquiries@aalaquis.com www.aalaquis.com

BARBADOS
Lot #4J, Lower Estate,
St Michael,
Barbados
Tel: 1 (246) 271 0180
Email: barbados@aalaquis.com

JAMAICA #48 Constant Spring Road, Kingston 10, Jamaica Tel: 1 (876) 648 3074 Email: jamaica@aalaquis.com