


## 2. Synthetic antibacterials

Data No. LA057-0000



Column	: Inertsil Ph 5µm 150 × 4.6 mm I.D.	Solutes	: 1. Olaquinox 2. Carbadox 3. Sulfadimidine 4. Thiamphenicol 5. Sulfadimidine 6. Sulfamonomethoxine 7. Furazolidone 8. Oxolinic acid 9. Sulfadimethoxine 10. Sulfaquinoxaline 11. Nalidixic acid 12. Piromidic acid 13. Nicarbazin
Column Cat. No.	: 5020-01301		
Eluent	: A: CH <sub>3</sub> CN B: 20 mM KH <sub>2</sub> PO <sub>4</sub> (pH 2.5; H <sub>3</sub> PO <sub>4</sub> ) A/B = 15/85 - 25 min - 50/50		
Flow rate	: 1.0 mL/min.		
Col. Temp.	: 40°C		
Detection	: UV 230 nm		
Injection Vol.	:		
Sample	: Synthetic Antibacterials		