


Analysis of Ascorbic acid [Vitamin C] (2)

Data No. LA494-0000



System	: GL-7400 system	Solutes	: 1. Erythorbic acid (5 mg/L)
Column	: Inertsil NH2 (5 μ m, 250 x 4.6 mm I.D.)		2. Ascorbic acid (5 mg/L)
Column Cat. No.	:		
Eluent	: A) CH ₃ CN B) CH ₃ OH C) 10 mM NaH ₂ PO ₄ D) 0.03 % Homocysteine A/B/C/D = 600/30/100/30 ,v/v/v/v		
Flow rate	: 1.0 mL/min		
Col. Temp.	: 40 °C		
Detection	: PDA 270 nm (GL-7452A PDA Detector)		
Injection Vol.	: 5 μ L		
Sample	:		