


Analysis of Glycyrrhizin

Data No. LA205-0000



Column : Inertsil ODS-P
5 μ m
150 \times 4.6mm I.D.

Solutes : 1. Glycyrrhizin (750ng)
2. p-Hydroxy benzoic acid propyl (150ng)

Column Cat. No. :

Eluent : H₂O / Acetic acid / MeCN (560 : 40 : 400)

Flow rate : 0.4 mL/min.

Col. Temp. : 25°C

Detection : UV 254 nm

Injection Vol. :

Sample : 1. Glycyrrhizin (750ng)
2. p-hydroxy benzoic acid propyl (150ng)