


## Analysis of Dansyl-Threonine (Thr)

Data No. LA384-0000



Column : Bioptic AV-1  
5 $\mu$ m  
150  $\times$  4.6 mm I.D.

Solutes : 1. Dansyl-D-Thr  
2. Dansyl-L-Thr

Column Cat. No. : 5020-02511

Eluent : Methanol / 50mM phosphate buffer (pH6.6) = 2 / 98 (w / w)

Flow rate : 1.0 mL/min.

Col. Temp. : 40°C

Detection : UV 254 nm

Injection Vol. : 2 $\mu$ L

Sample : Dansyl derivatives