


Analysis of Cyclodextrin

Data No. LA691-0000



Conditions

Column : Inertsil WP300 C18 (5 μm, 150 x 4.6 mm I.D.)
Column Cat. No. : 5020-05845
Eluent : A) CH₃CN
 : B) 25 mM phosphate buffer (pH 3)
 : A/B = 60/40, v/v
Flow Rate : 1.0 mL/min
Col. Temp. : 40 °C
Detection : RI
Injection Vol. : 20 μL

Analyte

1: α-Cyclodextrin
2: β-Cyclodextrin
3: γ-Cyclodextrin