


## Analysis of Benzene Derivatives

Data No. LA434-0000



Column : Inertsil Diol  
5  $\mu$ m  
150 x 4.6 mm I.D.

Column Cat. No. : 5020-05645

Eluent : A) CH<sub>3</sub>CN  
B) H<sub>2</sub>O  
A/B = 30/70, w/w

Flow rate : 1.0 mL/min.

Col. Temp. : 40°C

Detection : UV 254 nm

Injection Vol. :

Sample :

Solutes : 1. Acetophenone  
2. Benzene  
3. Toluene  
4. Naphthalene