


Analysis of Glycyrrhizin

Data No. La197-0000



Column : Inertsil Ph-3
5 μ m
150 \times 4.6mm I.D. Solutes : 1. Glycyrrhizin
2. n-Propyl paraben

Column Cat. No. :

Eluent : Acetonitrile / H₂O / Acetic Acid = 34 / 65 / 1

Flow rate : 1.0 mL/min.

Col. Temp. : 40°C

Detection : UV 254 nm

Injection Vol. :

Sample : Glycyrrhizin

p-Hydroxy benzoic acid n-Propyl ester