


## Organic solvents

Data No. GA079-0000



System : GC/FID  
Column : InertCap Pure-WAX  
0.25mm I.D. x 30m df=0.25µm  
Column Temp. : 40°C (5min) → 5°C/min → 200°C  
Carrier Gas : He 100kPa  
Injection : Split 1 : 80 0.1µl  
250°C  
Detection : FID Range 10°  
260°C  
Sample : Organic solvents

Solutes : 1. n-Propyl acetate  
2. Trichloroethylene  
3. iso-Butyl acetate  
4. Tetrachloroethylene  
5. iso-Butanol  
6. Methylcellosolve  
7. Ethylcellosolve  
8. Methylcellosolve acetate  
9. Ethylcellosolve acetate  
10. n-Butylcellosolve