


InertSearch™ for GC

InertCap® Applications

Isothiocyanate


Data No. GA114-0875



Conditions

System : GC/FID
Column : InertCap Pure-WAX
0.25 mm I.D. x 30 m df = 0.25 µm
Col. Cat. No. : 1010-68142
Col. Temp. : 40 °C (5 min hold) - 5 °C/min - 160 °C
Carrier Gas : He 100 kPa
Injection : Split flow 67 mL/min
200 °C
Detection : FID Range 10⁰
200 °C
Sample Size : 1000 µg/mL in Hexane
1 µL

Analyte : 1. Allyl isothiocyanate



Structures are created using Chemistry 4-D Draw which is provided by ChemInnovation Software, Inc.