


Hippuric acids

Data No. LA222-0000



Column : Inertsil ODS-2
5 μ m
150 \times 4.6mm I.D.

Column Cat. No. : 5020-01101

Eluent : CH₃CN/10mM KH₂PO₄+0.5% CH₃COOH
/ β -Cyclodextein=100/900/5

Flow rate : 1.0 mL/min.

Col. Temp. : 40°C

Detection : UV 227 nm

Injection Vol. :

Sample : Metabolite

Solutes : 1. Mandelic acid
2. Hippuric acid
3. o-Methyl hippuric acid
4. p-Methyl hippuric acid
5. m-Methyl hippuric acid
each 10 μ g