


Analysis of Benzene Derivatives

Data No. LA453-0000



Column : Inertsil WP300 Diol
5 μ m
250 x 4.6 mm I.D.

Column Cat. No. : 5020-05946

Eluent : A) CH₃CN
B) H₂O
A/B = 24.5/75.5, w/w

Flow rate : 1.0 mL/min.

Col. Temp. : 40°C

Detection : UV 254 nm

Injection Vol. :

Sample :

Solutes : 1. Acetophenone
2. Benzene
3. Toluene
4. Naphthalene