


Analysis of Fluocinonide

Data No. LA169-0000



Rs=12.2

Column : Inertsil ODS-3V
5 μ m
250 \times 4.6mm I.D.

Solutes : 1. Fluocinonide(0.2mg/mL)
2. n-Propylbenzoate(8 μ L/mL)

Column Cat. No. :

Eluent : CH₃CN/H₂O=50/50

Flow rate : 2.0 mL/min.

Col. Temp. : 40°C

Detection : UV 254 nm

Injection Vol. : 20 μ L

Sample : 1. Fluocinonide(0.2mg/mL)
2. n-Propylbenzoate(8 μ L/mL)